MR 36 Auto Hematology Analyzer

Operator's Manual

CE

Preface

Thank you for purchasing the Auto Hematology Analyzer manufactured by MR.

Read and understand the entire operator's manual before operating this device. Store this operator's manual properly for future reference.

Product name: Auto Hematology Analyzer

Model:MR36

Product Components: Blood Aspiration Module, Dilution Unit, Cleaning Unit, Analyzing and Measuring Unit and Microprocessor.

Scope of Use: blood cell counting, white blood cell 3-part classification and hemoglobin concentration measurement in clinical examinations.

Copyright

© MR International Healthcare Technology Co., Ltd. All rights reserved. This document contains proprietary information of MR International Healthcare Technology Co., Ltd. (hereinafter referred to as MR) No part of this document may be reproduced, copied, modified, disclosed, or transmitted in any form or by any means without prior written consent of MR,

This document is intended for users of MR equipment, who are authorized to use this document as they purchase MR equipment. Unauthorized persons are not allowed to use this document.

All information in this document does not constitute a warranty of any kind, express or implied, including but not limited to, the implied warranties of merchantability and fitness for a particular purpose. Every effort has been made in the preparation of this document to ensure accuracy of the contents. However, MR assumes no liability or responsibility for any errors or omissions in the contents of this document. MR reserves the right to improve any product at any time to enhance product reliability, functionality, or design.

Declaration

This operator's manual may be modified without notice.

MR reserves the right of final interpretation of this operator's manual.

The pictures in this operator's manual are for reference only. If there is inconsistency between the pictures and the actual product, the actual product shall prevail. Do not use the pictures for other than intended use.

MR shall be responsible for the safety, security, and performance of the product only when all of the following conditions are met:

- The assembly, re-commissioning, extension, modification, and repair of the product are performed by the authorized personnel of MR.
- The product is operated based on this operator's manual.
- The electrical appliances in the relevant working room comply with applicable national and local requirements.

Contact Info for After-Sales Services

MR International Healthcare Technology Co., Ltd.

Address:Room 1501,Grand Millennium Plaza,181 Queen's Road Central,HK

Tel:+8522 2158 1807 Fax:+852 3011 6829

E-mail:sales@mr-healthcare.com service@mr-healthcare.com

Website:http://www.mr-healthcare.com

Contents

Preface	i
1 Manual Overview	1
1.1 Introduction	1
1.2 Who Should Read This Manual	1
1.3 How to Find Information	1
1.4 Conventions Used in This Manual	2
1.5 Symbol Conventions	2
1.6 Safety Information	5
2 Installation	6
2.1 Introduction	6
2.2 Installation Personnel	6
2.3 Installation Requirements	6
2.4 Damage Inspection	7
2.5 Unpacking	8
2.6 Connecting the Analyzer System	8
2.6.1 Electrical Connections	8
2.6.2 Reagent Connections	9
2.6.3 Installing the Diluent Float Sensor and Replacing the Reagents	10
2.6.4 Installing the Waste Float Sensor	11
2.7 Installing Thermal Paper	12
3 System Overview	14
3.1 Introduction	14
3.2 Intended Use	14
3.3 Measurement Parameters	14
3.4 Structure of the Analyzer	15
3.4.1 Host	16
3.4.2 Touch Screen	18
3.4.3 Aspirate key	18
3.4.4 Power/Status indicator	18
3.4.5 Paper Feed Key	18
3.4.6 Thermal Printer	18
3.4.7 Power Switch	19
3.4.8 USB interface	19
3.4.9 Network interface	
3.4.10 External Equipment (Optional)	
3.5 User Interface	
3.6 Reagents, Controls and Calibrators	21

	3.6.1 Reagents	. 21
	3.6.2 Controls and Calibrators	. 22
4 ۱	Vorking Principle	. 23
	4.1 Introduction	. 23
	4.2 Aspiration	. 23
	4.3 Dilution	. 23
	4.3.1 Dilution Procedures in Whole-Blood Mode	23
	4.3.2 Dilution Procedures in Predilute Mode	24
	4.4 WBC/RBC/PLT Measurement	. 25
	4.4.1 Electrical Impedance Method	. 25
	4.4.2 Derivation of WBC-Related Parameters	. 26
	4.4.3 RBC	. 27
	4.4.4 PLT	. 28
	4.5 HGB Measurement	. 28
	4.5.1 Colorimetric Method	. 28
	4.5.2 HGB	. 28
	4.6 Flushing	. 29
5 \$	Setup	. 30
	5.1 Introduction	. 30
	5.2 Interface Introduction	. 30
	5.3 System Settings	. 31
	5.3.1 Date and Time	. 31
	5.3.2 Input Settings	. 32
	5.3.3 Lab Information	. 33
	5.3.4 Auto Maintenance	. 35
	5.4 Parameter Settings	. 35
	5.4.1 Data Dictionary	. 35
	5.4.2 Parameter Unit	. 39
	5.4.3 Ref. Range	. 41
	5.4.4 Microscopic Exam. Settings	. 48
	5.4.5 Customized Parameters	. 50
	5.5 Meterage Settings	. 53
	5.5.1 Gain Settings	. 53
	5.5.2 Flag	. 54
	5.6 Communication	
	5.6.1 Host Network Settings	
	5.6.2 LIS Communication	. 57
	5.7 User Management	
	5.7.1 Accessing the Interface	
	5.7.2 Creating a User	
	5.7.3 Editing a User	
	5.7.4 Deleting a User	
	5.7.5 Setting the Default User	
	5.7.6 Changing Password	
	5.7.7 Resetting Password	. 64

5.8 Print Settings	65
5.9 Auxiliary Settings	68
6 Daily Operations	71
6.1 Introduction	71
6.2 Pre-operation Preparation	71
6.3 Startup	72
6.4 Daily Quality Control	73
6.5 Sample Collection and Handling	74
6.5.1 Running the Venous Whole Blood Samples	74
6.5.2 Capillary Whole Blood Samples	75
6.5.3 Prediluted Samples	75
6.6 Sample Analysis	77
6.7 Shutdown	77
7 Sample Analysis	79
7.1 Introduction	79
7.2 Interface Introduction	79
7.3 Entering Sample Information	81
7.4 Running Samples	85
7.5 Dealing with the Analysis Results	86
7.5.1 Automatic saving of analysis results	86
7.5.2 Parameter Flags	87
7.5.3 Flags of Abnormal Blood Cell Differential or Morphology	87
7.6 Functions of the Buttons	88
7.6.1 Previous/Next	88
7.6.2 Mode & ID	88
7.6.3 Pre-entry	88
7.6.4 Validate/Cancel Validation	88
7.6.5 Print	89
7.6.6 Customized Parameters	89
7.6.7 Microscopic Exam. Parameters	89
7.6.8 Communication	
7.6.9 Edit Result	
7.6.10 Delete	
8 Result Review	94
8.1 Introduction	94
8.2 Interface Introduction	94
8.3 Sample List	95
8.4 Functions of the Buttons	
8.4.1 Validate	
8.4.2 Cancel Validation	
8.4.3 Print	
8.4.4 Delete	
8.4.5 Export	
8.4.6 Edit Result	
8.4.7 Patient Info.	102

8.4.8 Query	106
8.4.9 Graph	107
8.4.10 Microscopic Exam Para.	108
8.4.11 Customized Parameters	109
8.4.12 Run Chart	110
9 Quality Control	114
9.1 Introduction	114
9.2 L-J Quality Control	114
9.2.1 QC Principle	114
9.2.2 QC Settings	115
9.2.3 Quality Control Analysis	119
9.2.4 QC Result Review	125
9.3 X-B Quality Control	135
9.3.1 QC Principle	135
9.3.2 QC Settings	135
9.3.3 Quality Control Analysis	138
9.3.4 QC Result Review	138
10 Calibration	147
10.1 Introduction	147
10.2 When to Calibrate	147
10.3 How to Calibrate	147
10.3.1 Preparation	148
10.3.2 Manual Calibration	149
10.3.3 Auto Calibration Using Calibrators	151
10.3.4 Auto Calibration Using Fresh Blood Samples	153
10.4 Verifying Calibration Coefficients	154
11 Reagent Management	155
11.1 Accessing the Interface	155
11.2 Reagent Information Settings	156
11.2.1 Open system	156
11.2.2 Closed system	158
11.3 Reagent Replacement	160
12 Service	162
12.1 Introduction	162
12.2 Maintenance	162
12.2.1 Reagent Replacement	163
12.2.2 Cleaning	165
12.2.3 Maintenance	166
12.2.4 Comprehensive Device Maintenance	170
12.2.5 Auto Clean	175
12.2.6 Auto Prompt for Cleanser Soak	175
12.2.7 Auto Sleep	176
12.3 Self-inspection	176
12.3.1 Syringe and Sampling Mechanism	176

12.3.2 Pressure and Vacuum	177
12.3.3 Valve & Pump	178
12.3.4 Others	178
12.4 System Status	179
12.4.1 Temperature	179
12.4.2 Voltage and Current	180
12.4.3 Version Information	180
12.5 Log	182
12.5.1 All Logs	183
12.5.2 Parameter Revision Logs	183
12.5.3 Fault Logs	184
12.5.4 Other Logs	185
12.6 Touch Screen Calibration	186
13 Troubleshooting	187
13.1 Introduction	187
13.2 Dealing with Error Messages	187
13.3 Error Message Reference	188
Appendix A Specifications	192
A.1 Classification	192
A.2 Reagents	192
A.3 Parameters	192
A.4 Performance Specifications	193
A.4.1 Display Range	193
A.4.2 Normal Background	193
A.4.3 Linearity Range	194
A.4.4 Repeatability	194
A.4.5 Carryover	194
A.5 Input/output Device	195
A.6 EMC Description	195
A.7 Environmental conditions	195
A.8 Dimensions and Weight	196
A.9 Expected service life	
A.10 Contraindications	196
Appendix B Terms and Abbreviations	197
Annendix C Packing List	198

1 Manual Overview

1.1 Introduction

This chapter explains how to use this operator's manual of Auto Hematology Analyzer, which is shipped with the auto hematology analyzer and contains reference information about the analyzer and procedures for operating, troubleshooting and maintaining the analyzer.

Read this manual carefully before operating the analyzer and operate your analyzer in strict accordance with this manual.

1.2 Who Should Read This Manual

This manual contains information written for clinical laboratory professionals to:

- Learn about the hardware and software of the analyzer.
- Customize system settings.
- Perform daily operations.
- Perform system maintenance and troubleshooting.

1.3 How to Find Information

This operator's manual comprises 13 chapters and 3 appendix. Find the information you need by referring to the table below.

See	You can find
1 Manual Overview	Instructions for using the auto hematology analyzer.
2 Installation	Installation requirements for the auto hematology analyzer.
3 System Overview	Applications, measurable parameters, instrument configuration, software interface and software operations of the auto hematology analyzer.
4 Working Principle	Measuring principle and procedures of the auto hematology analyzer.
5 Setup	Settings of the system parameters such as the software date format and parameter units.
6 Daily Operations	Daily operations such as sample collection and preparation, the analysis procedures, startup and shutdown of the instrument.
7 Sample Analysis	Sample analysis procedure and handling of the analysis results.

See	You can find
8 Result Review	Review of the analysis results.
9 Quality Control	Basic requirements for quality control and the quality control methods provided by the auto hematology analyzer.
10 Calibration	Basic requirements for calibration and the calibration methods provided by the auto hematology analyzer.
11 Reagent Management	Settings and management of the reagents for the auto hematology analyzer.
12 Service	Methods for maintaining and testing the auto hematology analyzer.
13 Troubleshooting	Troubleshooting methods for the auto hematology analyzer.
Appendix A Specifications	Specification indicators of the auto hematology analyzer.
Appendix B Terms and Abbreviations	Terms and abbreviations of the auto hematology analyzer.
Appendix C Packing List	Packing list of the auto hematology analyzer.

1.4 Conventions Used in This Manual

The texts with special meaning in the Manual are highlighted by different fonts and formats.

Format	Definition
[XX]	All uppercase characters enclosed in [] indicate the name of a key on the analyzer or the peripheral keyboard, such as [ENTER].
XX	Bold characters indicate text displayed on the screen, such as Report .
XX	XX indicates variables and the specific content depends on the actual situation.
XX	Bold and italic characters Indicate chapter titles, such as 1.1 Introduction.

1.5 Symbol Conventions

The following symbols are used to indicate danger and alert messages in this manual.

When you see	Then
	Follow the instruction below the symbol to avoid potential biocontamination.

2 Co

When you see	Then
warning	Follow the instruction below the symbol to avoid personnel injury.
A CAUTION	Follow the instruction below the symbol to avoid analyzer damage and failure, or unreliable analysis results.
NOTE	Follow the instruction below the symbol. The symbol highlights the important information in operating procedures that calls for special attention.
<u> </u>	Puncture Warning: The sampling probe is sharp and may contain biohazardous materials. Special care should be taken when working with it.

The analyzer or the outer packaging may have the following labels or symbols.

NOTE

- If the labels are damaged or missing, please contact MR or MR's agents for replacement.
- All illustrations in this manual are provided as references only. They may not necessarily reflect actual analyzer configuration or display.

When you see	It means
<u> </u>	Caution
	Biohazard
	Exercise caution to prevent puncture
	Instruction for Moving
	Network interface
	Protective grounding
\sim	Alternating current (AC)
IVD	For in vitro diagnosis only

When you see	It means
LOT	Lot No.
	Expiry date
SN	Serial No.
C€	The device is in full compliance with the council directive concerning in vitro diagnostic medical devices 98/79/EC.
EC REP	Authorized Representative in the European Community
	Date of manufacture
	Manufacturer
40°C	Storage temperature
100	Humidity level for storage
KPa 106	Atmospheric pressure level for storage
\bigcap i	Consult the operator's manual
**	Avoid sunlight
Ť	Keep dry
	No rolling
	No Stacking.
<u> </u>	Let this side face upward.

When you see	It means
1	Fragile, handle with care
	Recyclable materials
X	The analyzer, after being scrapped, should not be disposed with other household garbage, instead, it should be collected and recycled following the disposal instructions for scrapped electronic and electrical equipment.

1.6 Safety Information



- All the samples, controls, calibrators, reagents, wastes and areas in contact with them are
 potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms,
 etc.) and follow laboratory safety procedures when handling relevant items and areas in the
 laboratory.
- If leak happens to the analyzer, the leak liquid is potentially biohazardous.



WARNING

- Please check the firmness of all the door/ covers/panels before running the analyzer to prevent unexpected opening or loosening when the analyzer is working.
- Make sure all the safety measures are taken. Do not disable any safety device or sensor.
- Please respond to any alarm and error message immediately.
- Do not touch the moving parts.
- Contact MR or MR authorized agents upon the identification of any damaged part.
- Be careful when opening/closing and removing/installing the doors, covers and panels of the analyzer.
- Dispose the analyzer according to government regulations.



CAUTION

- Please use the analyzer in strict accordance with this manual.
- Please take proper measures to prevent the reagents from being polluted.

2 Installation

2.1 Introduction



WARNING

Installation by personnel not authorized or trained byMR International Healthcare Technology Co.,Ltd maycause personal injuryor damage to the analyzer. Do not install the analyzer without the presence of MR- authorized personnel.

Your analyzer has passed strict tests before it is shipped from the factory. Internationally-recognized symbols and instructions show the carrier how to properly handle this electronic instrument in transportation. When you receive your analyzer, carefully inspect the packaging. If you see any sign of mishandling or damage, contact MR customer service department or your local agent immediately.

2.2 Installation Personnel

The analyzer should only be installed by MR or its authorized agents. You need to provide the appropriate environment and space. When the analyzer needs to be relocated, please contact MR or your local agents.

When you receive the analyzer, please notify MR or your local agent immediately.

2.3 Installation Requirements



WARNING

- Connect only to a properly grounded outlet.
- Before turning on the analyzer, make sure the input voltage meets the requirements.



CAUTION

- Using a patch board may introduce electrical interference and generate incorrect analysis results.
 Please place the analyzer near the electrical outlet to avoid using the patch board.
- Please use the original electrical wires shipped with the analyzer. Using other electrical wires may damage the analyzer or generate incorrect analysis results.

Installation requirements for the analyzer are as follows.

Installation Environment	Requirements	
	Level ground and stable workbench with load capacity ≥50kg.	
	Free of dust, mechanical vibration, heat and wind sources, contamination, heavy-noise source or electrical interference.	
Site	Avoid direct sunlight and keep good ventilation.	
Sito Sito	It's recommended to evaluate the electromagnetic environment of the laboratory before operating the analyzer.	
	Keep the analyzer away from sources of strong electromagnetic interference, otherwise, its proper functioning may be affected.	
	 At least 50cm from each side, which is the preferred access to perform service procedures. 	
Space (In addition	At least 20cm from the back for cabling and ventilation.	
to the space required for the analyzer itself, set aside)	 Enough room on and below the countertop to accommodate for the diluent and waste containers. 	
	 Place the analyzer near the electrical outlet and avoid being blocked by any objects, so that you can disconnect the power plug easily as required. 	
Optimal operating temperature	10°C~30°C	
Optimal operating humidity	20%~85%	
Operating atmospheric pressure	70kPa~106kPa	
Ventilation	Keep air exchange to ensure good air circulation. The wind should not blow directly at the analyzer.	
Power	AC100V~240V, Input Power ≤200VA, 50/60HZ.	
Electromagnetic Wave	Keep the analyzer away from electric-brush motors, flashing fluorescent and electric-contact equipment which is switched on/off frequently.	
Waste Disposal	Dispose of the waste as per the requirements of the local environment protection authorities.	

2.4 Damage Inspection

Before packing and shipping, MR has applied rigid inspection on all the analyzers. Upon receiving the analyzer, please check carefully before unpacking to see if there are any of the following damages:

- The outer packaging is placed upside down or distorted.
- The outer packaging shows obvious signs of having been exposed to humid conditions.
- The outer packaging shows obvious signs of having been crashed.
- The outer packaging shows signs of having been opened.

Once you find the above damages, please notify your local agent immediately.

If the packaging is intact, please open the packaging in the presence of personnel from MR or its agents and apply the following inspections:

- Check if all the items listed in the packing list are in the packaging.
- Carefully inspect the appearance of all the items to check if they are damaged or distorted.

2.5 Unpacking

Please unpack the analyzer by taking the following steps:

- 1. Open the outer packing box; take out the accessory pack; take out the analyzer together with the protective and cushioning materials.
- 2. Remove the foam and the protective PE bag.
- Open the right door (open the linear-shaped cam lock on the right door with a slotted screwdriver).
- 4. Remove the binder clips, which are used for fixating two conveyor belts.
 - To avoid the possible collision resulting from the slippage caused by shaking and slanting during transportation, the central position of those two belts is fixated with binder clips before they are shipped from the factory. The binder clips must be removed during unpacking.
- 5. Remove the binder clips, which are used for fixating sampling assembly.

NOTE


To avoid damage during the transportation, the sampling assembly of the analyzer is fixated with clamps. Do remove the clamps before using the analyzer.

2.6 Connecting the Analyzer System

2.6.1 Electrical Connections

Please refer to Figure 2-1 for the electrical connections of the analyzer.

Figure 2-1 Connecting the electric lines



2.6.2 Reagent Connections



WARNING

- Be sure to dispose of reagents, waste, samples, consumables, etc. according to local legislations and regulations.
- The reagents are irritating to eyes, skin and mucosa. Wear proper personal protective
 equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when
 handling them in the laboratory.
- If the reagents accidentally spill on the skin, wash them off with plenty of water and if necessary, go see a doctor; if the reagents accidentally spill into the eyes, wash them off with plenty of water and immediately go see a doctor.




CAUTION

- Please make sure the length of the diluent pipe and the waste pipe should be no longer than 1500mm; the length of the lyse pipe and the cleanser pipe should be no longer than 850mm.
- Tighten the panel connector of the fluidic line so that the overall fluidic line is closed to prevent leakage and seepage caused by siphonage, etc.

Refer to Figure 2-2 for the connection of the reagents placed outside the analyzer.

Figure 2-2 Connecting reagents placed outside the analyzer



Refer to Figure 2-3 for the connection of the reagent placed inside the analyzer.

LYSE

Figure 2-3 Connecting reagents placed inside the analyzer (left door opened)

2.6.3 Installing the Diluent Float Sensor and Replacing the Reagents


Please install the diluent float sensor and replace the diluent as per the approaches stated in this section.

2.6.3.1 Installing the Diluent Float Sensor

Install the diluent float sensor according to the following steps.

- 1. Press down and remove the round cardboard with dotted cutting line on the top side of the diluent box so as to reveal a round hole.
- 2. Pull out the cover of the container so that the cardboard around the round hole can seize the neck under the vial cap to prevent invagination.
- 3. Turn and open the cap (keep the cap) and prevent any foreign objects from getting into the container.
- 4. Install the diluent float sensor assembly in the accessory pack as shown in Figure 2-4. The float sensor shall be kept as vertical as possible during installation and the self-contained cap of the sensor shall be tightened.

Figure 2-4 Installing the Diluent Float Sensor



2.6.3.2 Replacing Reagents

Steps for the replacing the diluent are the same as that for installing the sensor. Please keep the empty diluent container and the cap for future use.


2.6.4 Installing the Waste Float Sensor

NOTE

The float sensors used in the analyzer are only applicable to MR- supplied waste containers or the containers with the same specification and model (such as the vacant diluent container).

- 1. Take a proper waste container (it can be a vacant diluent container, the opening of which is required to be pulled out of the hole of the box to expose the opening) and open the vial cap.
- 2. Install the waste float sensor assembly in the accessory pack as shown in Figure 2-5. The float sensor shall be kept as vertical as possible during installation and the self-contained cap of the sensor shall be tightened at the same time to prevent the spilling of the waste.

Figure 2-5 Installing the Waste Float Sensor



The waste container can be replaced according to the steps mentioned above. The replaced waste shall be properly disposed to avoid contamination.



WARNING

Be sure to dispose of reagents, waste, samples, consumables, etc. according to local legislations and regulations.

2.7 Installing Thermal Paper



CAUTION

- Use only specified thermal paper. Otherwise, it may cause damage to the thermal printer head, or the printer may be unable to print, or poor print quality may result.
- Never pull the thermal printer paper with force when a recording is in process. Otherwise, it may cause damage to the thermal printer.
- Do not leave the thermal printer door open unless you are installing paper or removing error.
- Improper installation of thermal printer paper may jam the paper and/or result in blank printout.

NOTE

Remove the protective paper between the thermal printer head and the roller inside the thermal printer before installing thermal paper for the first time.


Follow the procedure below to install the thermal paper.

1. Use the latch (as shown in Figure 2-6) at the upper right corner of the thermal printer door to pull the door open.



Figure 2-6 Installing Thermal Paper (1)

2. Insert a new roll into the compartment as shown below.





3. Close the thermal printer door.

4. Check if paper is installed correctly and the paper end is feeding from the top.

Figure 2-8 Installing Thermal Paper (3)



5. To ensure the normal use of the thermal paper, press the feed key to start paper feeding, and then press the feed button again to stop feeding when a short paper is sent out.

3 System Overview

3.1 Introduction

Auto Hematology Analyzer is a quantitative, automated hematology analyzer and 3-part differential counter used in clinical laboratories.

This section describes in details the intended use, measurement parameters, structure, user interface and compatible reagents of the analyzer.

3.2 Intended Use

It's intended for blood cell counting, 3-part classification of white blood cell and hemoglobin concentration measurement in clinical examinations.

NOTE

The analyzer is intended for screening in the clinical examination. When making clinical judgment based on the analysis results, the doctors should also take into consideration the clinical examination results or other test results.

3.3 Measurement Parameters

As shown below, the analyzer provides quantitative analysis results for 21 hematology parameters and three histograms.

Туре	Parameter Name	Abbreviation
WBC	White Blood Cell count	WBC
(7 items)	Percentage of Granulocytes	Gran%
	Percentage of Lymphocytes	Lym%
	Percentage of Mid-sized Cells	Mid%
	Number of Granulocytes	Gran#
	Number of Lymphocytes	Lym#
	Number of Mid-sized Cells	Mid#
RBC (8 items)	Red Blood Cell count	RBC
	Hemoglobin Concentration	HGB

Туре	Parameter Name	Abbreviation
	Mean Corpuscular Volume	MCV
	Mean Corpuscular Hemoglobin	MCH
	Mean Corpuscular Hemoglobin Concentration	MCHC
	Red Blood Cell Distribution Width - Coefficient of Variation	RDW-CV
	Red Blood Cell Distribution Width - Standard Deviation	RDW-SD
	Hematocrit	нст
PLT	Platelet count	PLT
(6 items)	Mean Platelet Volume	MPV
	Platelet Distribution Width	PDW
	Plateletcrit	PCT
	Platelet-large cell ratio	P-LCR
	Platelet-large cell count	P-LCC
Histogram	White Blood Cell Histogram	WBC Histogram
	Red Blood Cell Histogram	RBC Histogram
	Platelet Histogram	PLT Histogram

3.4 Structure of the Analyzer



WARNING

- Please check the firmness of all the doors, covers and boards before running the analyzer.
- The analyzer is heavy, so moving by one person alone may cause injury. It is advisable for two
 people to move it together when the transportation is necessary, and make sure you follow the
 instructions and use the proper tools.
- Connect only to a properly grounded outlet.
- To avoid electrical shocks, disconnect the power supply before opening the cover.
- To prevent fire, use the fuses with specified model number and working current.




The sampling probe is sharp and may contain biohazardous materials. Special care should be taken when working with it.

3.4.1 Host

The Auto Hematology Analyzer consists of the main unit (analyzer) and accessories. The main unit is the main part for analysis and data processing.


Front of the analyzer

Figure 3-1 Front of the analyzer



- 1: Touch screen
- 3: Sample probe
- 5: Printer status indicator
- 2: Power/Status indicator
 - 4: Aspirate key
 - 6: Paper feed key
- 7: Key for opening the paper compartment of the thermal printer
- 8: Paper compartment of the thermal printer
- Back of the analyzer

Figure 3-2 Back of the analyzer




1: Power switch

- 2: AC input
- 3: Diluent presence detection connector
- 4: Waste level detection connector

5: DIL-E diluent inlet

- 6: Waste outlet
- Right side of the analyzer (right door opened)

Figure 3-3 Right side of the analyzer




- 1: Power Supply
- 3. Count bath assembly
- 5: USB interface
- 7: Pressure chamber
- 2: Sampling assembly
 - 4: Valve assembly
- 6: Network interface
 - 8: Pump
- Left side of the analyzer (left door opened)



WARNING

To prevent injuries, do not place your hands near the bottom guide tracks of the syringes when the analyzer is running.

Figure 3-4 Left side of the analyzer



- 1: Valve assembly
- 2: Liquid level detection unit

3: Syringe

3.4.2 Touch Screen

The touch screen is located on the front side of the analyzer for performing interface operations and displaying the information.

3.4.3 Aspirate key

The aspirate key is located in the middle of the front side (behind the sample probe) to start the sample analysis, to add diluent, or to cancel sleep.

3.4.4 Power/Status indicator

The status indicator is located in the middle section of the right part of the analyzer (front side). It shows the status of the analyzer including ready, running, error, sleep and on/off, etc.

The indicators change with the status of the main unit. Details are shown in Table 3-1.

Instrument Status Indicator Status Remarks Shutdown Off The main unit has been shut down. Stopped running with Red light on Stopped running with the occurrence of errors error conditions Running with error Red light flickering Running with the occurrence of errors conditions Time sequence Yellow light on Initialization or sleep status irrelevant to deactivated running Running Green light flickering Execution of the sequence actions is in process. Ready Green light on Execution of the sequence actions is allowed.

Table 3-1 Main Unit Status Indicators

NOTE

While the analyzer is running, if the indicator turns dim or off, please contact MR or MR's agent for maintenance.

3.4.5 Paper Feed Key

The paper feed key is located below the touch screen. After you press it, the built-in thermal printer will send out the paper with records.

3.4.6 Thermal Printer

The thermal printer is located below the touch screen. It will send out the paper with records after you press the paper feed key.

3.4.7 Power Switch



CAUTION

To avoid damage, do not power on/off the analyzer repetitively within a short time.

A power switch is located in the bottom back of the analyzer. It turns on or shuts down the analyzer.

3.4.8 USB interface

The USB interface is located on the right side of the main unit. There are 4 interfaces in total for external equipment (printer, barcode scanner, mouse or keyboard, and so on) connection or data transmission.

3.4.9 Network interface

The network interface is located on the right side of the main unit. There is 1 network interface in total for connecting with the Ethernet.

3.4.10 External Equipment (Optional)

The analyzer can be connected with the following external equipment:

Keyboard

The keyboard is connected with the USB interface on the right side of the analyzer for controlling the analyzer.

Mouse

The mouse is connected with the USB interface on the right side of the analyzer for operations on the analyzer.

Printer

The printer is connected with the USB interface on the right side of the analyzer for printing reports and other information displayed on the screen.

Barcode Scanner

The barcode scanner is connected with the USB interface on the right side of the analyzer for entering barcode information in an easy and fast way.

USB flash disk

The USB flash disk is connected with the USB interface on the right side of the analyzer for exporting sample data.

3.5 User Interface

After the startup procedure, you will enter the user interface (**Sample Analysis** as default). See Figure 3-5.



Figure 3-5 User Interface

The interface can be divided into several areas as follows according to their functions:

1 - Menu navigation area

On the top of the screen is the menu navigation area. Once a menu button is pressed, the system goes immediately to the corresponding screen.

2 - Menu content display area

It displays the selected screen and the corresponding function buttons.

3 - Error message area

Upon the occurrence of a system failure, the corresponding error message will appear in this area. When there is more than one failure, the error message for the latest failure will appear in this area.

Click in this area, you can deal with the failures in the popup dialog box of troubleshooting help. For more information, see *13 Troubleshooting*.

4 - Status display area

On the top right of the screen is the status display area where the connection status between the computer and the LIS system and printer status are displayed from left to right. The icons change with the status of the main unit, as shown in Table 3-2.

Table 3-2 Status Icon Description

Status	Icon	Remarks
LIS/HIS status	Gray icon	The computer is not connected to the LIS/HIS.

Status	Icon	Remarks
	Black icon	The computer is connected to the LIS/HIS.
Print status	Gray icon	The printer is not connected to the analyzer yet.
	Color icon	The printer is connected to the analyzer.

5 - Information area of the next sample

This area displays the information about the sample ID, sample position and blood mode of the next sample.

- 6 Username of the current user
- 7 Current date and time of the analyzer.

3.6 Reagents, Controls and Calibrators

Because the analyzer, reagents, controls, and calibrators are components of the system, system performance depends on the combined integrity of all the components. You should only use the MR specified reagents (see *A.2 Reagents*), which are formulated specifically for the fluidic system of your analyzer in order to achieve optimal system performance. Do not operate the analyzer using reagents from multiple suppliers. Under such circumstances, the analyzer may not achieve the performance specified in this manual and may generate unreliable results. All references to "reagents" in this manual refer to the reagents specifically formulated for this analyzer.

Each reagent package should be examined before use. Inspect the package for signs of leakage or moisture. Product integrity may be compromised in packages that have been damaged. If there is evidence of leakage or improper handling, do not use the reagent.

NOTE

- After long-distance transportation, the reagent must be allowed to settle for more than one day before use.
- Store and use the reagents by following the instructions for use of the reagents.
- When you have changed the diluents or lyses, run a background check to see if the results meet the requirement.
- Pay attention to the expiration dates and open-container stability days of all the reagents. Be sure not to use expired reagents.

3.6.1 Reagents

The following reagents are intended to be used with the analyzer for 3-part diff counting, daily cleaning and other operations.

DIL-E Diluent

This product is intended for sample dilution and preparation of cell suspension before running the samples.

LYE-1 LYSE

This product is intended for lysing the red blood cells, determining the hemoglobin, white blood cell classification and counting the total number of white blood cells.

CLE-P Cleanser

This product is intended for cleaning the fluidic system of the analyzer and regular instrument cleaning.

3.6.2 Controls and Calibrators

The controls and calibrators are used for quality control and analyzer calibration.

The controls are commercially prepared whole-blood products used to verify that the analyzer is functioning properly. They are available in low, normal, and high levels. Daily use of all levels verifies the normal operation of the analyzer and ensures the acquisition of reliable results. The calibrators are commercially prepared whole-blood products used to calibrate the analyzer.

Read and follow the instructions to use the controls and calibrators.

The "calibrators" and "controls" mentioned in this manual refer to MR - specified calibrators and controls and need to be purchased from MR or its specified agent.

4 Working Principle

4.1 Introduction

The measurement methods used in this analyzer are: the Electrical Impedance method for determining the WBC, RBC and PLT and their volume distribution; the colorimetric method for determining the HGB. During each analysis cycle, the sample is aspirated, diluted and mixed before the determination for each parameter is performed.

4.2 Aspiration

The analyzer supports Whole Blood mode (including Venous Whole Blood and Capillary Whole Blood) and Predilute mode.

In Whole Blood mode, the analyzer will aspirate quantitative whole blood sample.

In Predilute mode, the analyzer will aspirate the prediluted sample (with the dilution ratio of 1:10) which is a mixture of 20μ L of whole blood/capillary blood sample and 180μ L of diluent the diluted sample thus prepared is then delivered to the analyzer for sampling and aspiration.

4.3 Dilution

After being aspirated into the analyzer, the sample is divided into two parts. After the reaction with reagents in parallel dilution procedures, each part forms the sample for red blood cell/platelet, white blood cell count/hemoglobin measurement.

To meet different needs, the analyzer offers two working modes –Whole Blood and Predilute. The dilution procedures for whole blood samples and prediluted samples will be presented on the following pages.

4.3.1 Dilution Procedures in Whole-Blood Mode

Dilution Procedures in Whole-Blood Mode are shown in Figure 4-1.




Figure 4-1 Dilution Procedures in Whole-Blood Mode

is the dilution procedure for red blood cell and platelet; is the dilution procedure for white blood cell count/hemoglobin; namely CBC.

4.3.2 Dilution Procedures in Predilute Mode

The dilution procedure for the prediluted sample is shown in Figure 4-2.




Figure 4-2 Dilution Procedures in Predilute Mode

① is the dilution procedure for red blood cell and platelet; ② is the dilution procedure for white blood cell count/hemoglobin; namely CBC.

4.4 WBC/RBC/PLT Measurement

The analyzer detects the white blood cell count, red blood cell count and platelet count and their volume distribution by impedance method and eventually obtains the results of related parameters.

4.4.1 Electrical Impedance Method

WBCs/RBCs/PLTs are counted and sized by the Electrical Impedance method. This method is based on the measurement of changes in electrical resistance produced by a particle, which in this case is a blood cell, suspended in a conductive diluent as it passes through an aperture of known dimensions. An electrode is submerged in the liquid on both sides of the aperture to create an electrical pathway. As each particle passes through the aperture, a transitory change in the resistance between the electrodes is produced. This change produces a measurable electrical pulse. The number of pulses thus generated is equal to the number of particles that passed through the aperture. The amplitude of each pulse is proportional to the volume of each particle.




Figure 4-3 Electrical Impedance method

Each pulse is amplified and compared to the internal reference voltage channel, which only accepts the pulses of a certain amplitude. If the pulse generated is above the WBC/RBC/PLT lower threshold value, it is counted as a WBC/RBC/PLT. The cell volume distribution is determined by the cell count within each channel classified by the pulse amplitude.

The analyzer presents the WBC/RBC/PLT histogram, where the x-coordinate represents the cell volume (fL) and the y-coordinate represents the number of the cells.

4.4.2 Derivation of WBC-Related Parameters

White blood cells have a variety of types and can be categorized according to their volume. The volume of each type of cells varies with the added diluent, lyse and the lysing time. With the action of reagents, white blood cells can be classified into three groups, in the order from small volume to large volume: Lymphocytes, Mid-sized Cells (including Monocytes, Eosinophils and Basophils) and Granulocyte.

Based on the white blood cell histogram and the analysis for the Lym zone, Mid zone and Gran zone, the analyzer can get the percentage of lymphocytes (Lym%), the percentage of mid-sized cells (Mid%) and the percentage of granulocytes (Gran%), and then get the number of lymphocytes (Lym#), the number of mid-sized cells (Mid#) and the number of granulocytes (Gran#) based on the calculation with the white blood cell count obtained with the electrical impedance method. The unit of the number of cells is 10^9 /L.

White Blood Cell count

WBC count is the number of leukocytes measured directly by counting the leukocytes passing through the aperture.

Percentage of Lymphocytes (Lym%)

 $Lym\% = \frac{Particle\ countin\ the\ Lym\ zone}{Sum\ of\ particle\ countin\ the\ Lym\ zone, Mid\ zone\ and\ Granzone} \times 100\%$

Percentage of Mid-sized Cells

Mid% = $\frac{\text{Particle countin the } \textit{Mid} \, \text{zone}}{\text{Sum of particle countin the } \textit{Lym} \, \text{zone}, \textit{Mid} \, \text{zoneand } \textit{Gran} \text{zone}} \times 100\%$

Percentage of Granulocytes (Gran%)

Gran% = $\frac{\text{Particle count in the } Granzone}{\text{Sum of particle count in the } Lym \text{ zone,} Mid \text{ zone and } Granzone} \times 100\%$

Number of lymphocytes (Lym#)

 $Lym# = WBC \times Lym\%$

Number of Mid-sized Cells

 $Mid# = WBC \times Mid\%$

Number of Granulocytes (Gran#)

 $Gran# = WBC \times Gran\%$

4.4.3 RBC

Red Blood Cell count

RBC (12¹²/L) is the number of erythrocytes measured directly by counting the erythrocytes passing through the aperture.

Mean Corpuscular Volume

Based on the RBC histogram, this analyzer calculates the mean corpuscular volume (MCV) and expresses the result in fL.

 Hematocrit (HCT), Mean Corpuscular Hemoglobin (MCH), Mean Corpuscular Hemoglobin Concentration (MCHC)

This analyzer calculates the HCT (%), MCH (pg) and MCHC (g/L) as follows, where the RBC is expressed in 10^{12} /L, MCV in fL and HGB in g/L.

$$HCT = \frac{RBC \times MCV}{10}$$

$$MCH = \frac{HGB}{RBC}$$

$$MCHC = \frac{HGB}{HCT} \times 100$$

• Red Blood Cell Distribution Width Coefficient of Variation (RDW-CV)

Based on the RBC histogram, this analyzer calculates the CV (Coefficient of Variation, %) of the erythrocyte distribution width.

• Red Blood Cell Distribution Width Standard Deviation (RDW-SD)

RDW-SD (RBC Distribution Width - Standard Deviation, fL) is obtained by calculating the

standard deviation of the red blood cell size distribution.

4.4.4 PLT

Platelet count (PLT count, 10⁹/L)

PLT is measured directly by counting the platelets passing through the aperture.

Mean Platelet Volume (MPV, fL)

Based on the PLT histogram, this analyzer calculates the MPV.

Platelet Distribution Width (PDW)

PDW is the geometric standard deviation (GSD) of the platelet size distribution. Each PDW result is derived from the platelet histogram data and is reported as 10(GSD).

Plateletcrit (PCT)

This analyzer calculates the PCT as follows and expresses it in %, where the PLT is expressed in 10⁹/L and the MPV in fL.

$$PCT = \frac{PLT \times MPV}{10000}$$

Platelet-Large Cell Count (P-LCC, 10⁹/L)

P-LCC is measured directly by counting the large platelets passing through the aperture.

Platelet-Large Cell Ratio (P-LCR)

$$P-LCR = \frac{P-LCC}{PLT} \times 100\%$$

4.5 HGB Measurement

HGB is determined by the colorimetric method.

4.5.1 Colorimetric Method

The WBC/HGB diluent is delivered to the HGB bath where it is mixed with a certain amount of lyse, which converts hemoglobin to a hemoglobin complex that is measurable at 525 nm. An LED is mounted on one side of the bath and emits a beam of monochromatic light with a central wavelength of 525nm. The light passes through the sample and is then measured by an optical sensor mounted on the opposite side. The signal is then amplified and the voltage is measured and compared with the blank reference reading (readings taken when there is only diluent in the bath).

4.5.2 HGB

The HGB is calculated using the following equation and expressed in g/L.

$$HGB(g/L) = Constant \times Ln \left(\frac{Blank\ Photocurrent}{Sample\ Photocurrent} \right)$$

4.6 Flushing

After each analysis cycle, each component of the analyzer is flushed.

5 Setup

5.1 Introduction

The analyzer has been initialized before delivery. The interfaces upon the initial startup of the analyzer are system settings by default. Some parameters of the analyzer can be reset to meet various demands in practical applications.

The analyzer divides the operators into two access levels, common user and administrator. Note that an administrator can access all the functions accessible to a common user. This chapter introduces how to customize your analyzer as an administrator.

5.2 Interface Introduction

After logging in the software system (see **6.3 Startup**), click and choose Setup to access the Setup interface.

See Figure 5-1.



Figure 5-1 Setup

The administrator is allowed to set the following functions in the **Setup** interface:

- System settings
- Parameter settings
- LIS communication

- Meterage settings
- User management
- Print settings
- Auxiliary settings

5.3 System Settings

5.3.1 Date and Time

You can set the current date and time, as well as the date display format in the analyzer system. Specific steps are shown below:

1. Click **Date and Time** in the **System** area.

The date and time format setting interface pops up.



2. Click the **Date and Time** dropdown list and set the current date and time of the system in the popup dialog box.




Related descriptions:

> The input sequence of the controls is the same with the date fomat on the top right corner of

the dialog box. For example, if the data format is **yyyy/MM/dd HH:mm:ss**, you should input the data in the sepuence of year, month, date, hour, minute, and second.

- > Click or to select the date and time or click the textbox to enter them directly.
- > Click to clear the data and input again.
- 3. Click **OK** to save and close the message box.
- 4. Select the format setting from the dropdown list of the **Date Format**. See Figure 5-2.


Figure 5-2 Setting the Date Format



5. Click Apply.

The system message will pop up, indicating the successful setting. See Figure 5-3.

Figure 5-3 Successful Setting of the Date Format



The date and time at the bottom right corner will be displayed in the newly set format as shown in 06/26/2015 15:11:44.

- 6. Click **OK** to close the message box.
- 7. Click **OK** to exit.

5.3.2 Input Settings

Click **Input setting** in the **System** area, and then you can set the soft keyboard for screen input. As shown in Figure 5-4, You can set to turn the soft keyboard on or off.

Figure 5-4 Input Settings




Soft Keyboard

Open (default)

You can enter content using the soft keyboard popped up on the screen. Functions and applications for the keys are shown in Figure 5-5.

Figure 5-5 Soft Keyboard



Close

You need to use an externally connected USB keyboard for entering content.

5.3.3 Lab Information

Click **Lab Information** in the **System** selection, then you can set the lab information. See Figure 5-6.




Figure 5-6 Setting Lab Information

NOTE

Only the administrator has the access for setting the lab information. General users are only allowed to browse such information.

Refer to the table below for the detailed instructions of parameter setting.

Table 5-1 Setting Lab Information

Parameter	Setting Description	
Hospital Name	Enter the name of the hospital where the lab is located.	
Lab Name	Enter the lab name.	
Responsible Person	Enter the responsible person of the lab.	
Contact Information	Enter the contact information (telephone number or E-Mail) of the lab.	
Contact in Service Department	Enter the name of the contact person in Service Department.	
Contact Information of Service Department	Enter the contact information of the contact person in the Service Department.	

Parameter	Setting Description	
Analyzer SN	Display the serial number of the analyzer. Read only.	
Installation Date	Display the installation date of the analyzer. Read only.	
Remarks	Enter the remarks regarding the lab.	

5.3.4 Auto Maintenance

Click **Auto Maintenance** in the **System** selection to access the **Auto Maintenance** setting interface. The system auto sleep waiting time and cleanser maintenance time can be set in the **Auto Maintenance** interface.

Auto Maintenance

Auto Sleep

Wait

60 minutes [15 120]

Auto Cleanser Soak

Start Time

17:00 minutes [0:00 23:59]

Figure 5-7 Auto Maintenance

Auto Sleep

In the **Wait** textbox, the administrator is allowed to set the waiting time for entering the sleep state after the main unit is halted. The range is between 15 and 120 minutes and the default value is 60 minutes.

OK

Cancel

Auto Cleanser Soak

The administrator is allowed to set the start time of the cleanser soak in the **Start Time** textbox. The acceptable value ranges from 0:00 to 23:59 and the default value is 17:00

5.4 Parameter Settings

5.4.1 Data Dictionary

You can set shortcut codes for the relevant items of the patient information.

Apply

If a shortcut code is set, the shortcut code corresponding to the above mentioned item can be entered directly when the information is input or numbered, then the complete information can be displayed without entering (or selecting) complete information.

It is a shortcut operation. Different items can share one shortcut code.

5.4.1.1 Accessing the Interface

Click **Data** in the **Para.** selection to access the data dictionary setting interface. See Figure 5-8. You can set the shortcut code for the relevant items of the patient information in this interface.




Figure 5-8 Shortcut Code

You can set the shortcut code for the following items: **Department**, **Submitter**, **Patient Type**, **Gender**, **Area**, **Bed No**. and **Sample Type**.

5.4.1.2 Adding a New Item


This section takes the adding of a new department as an example to introduce the method for adding a new item and its shortcut code. The method for adding other new items is similar and is not introduced in details herein.

Steps for adding a new department are shown as follows:

1. Click New in the Department interface.

A dialog box will pop up as shown in Figure 5-9.

Figure 5-9 Adding a New Item



2. Enter a new department name, shortcut code and remarks.

NOTE

- Newly added department name must be entered and it can not be the same as existing ones.
- The shortcut code is not necessary to be entered, but once set, every code must be unique.
- Click **OK** to save the information about the new department.
 Information about the newly added department will be displayed in the department interface. See Figure 5-10.

Data Name **Shortcut Code** Remarks Department 7 Internal Medicine Nk Submitter Wk Surgery OP Ophthalmology Patient Type A Gender Area Bed No. Sample Type New Edit Delete Cancel

Figure 5-10 Information of the Newly Added Department

5.4.1.3 Editing Items/Shortcut Code


This section takes the editing of a department as an example to introduce the method for editing items and its shortcut code. The method for editing other new items is similar and is not introduced in details herein.

Steps for editing a department are shown as follows:

1. Select the department to be modified in the **Department** interface (for example the Internal Medicine), then click **Edit**.

A dialog box will pop up as shown in Figure 5-11.

Figure 5-11 Editing Item/Shortcut Code



2. Modify the **Name**, **Shortcut Code** and **Remarks** in each textbox according to the actual demand.

NOTE

- Newly added department name must be entered and it can not be the same as existing ones.
- The shortcut code is not necessary to be entered, but once set, every code must be unique.
- Click OK to save the information.


5.4.1.4 Deleting a Shortcut Code

This section takes the deleting of a department as an example to introduce the method for deleting items and this shortcut code. The method for deleting other new items is similar and is not introduced in details herein.

Steps for deleting a department are shown as follows:

Select the department to be deleted in the **Department** interface, and then click **Delete**.
 A dialog box will pop up as shown below.

Figure 5-12 Deleting a Department



2. Click **Yes** to delete the department.

5.4.2 Parameter Unit

Some of the parameters of the analyzer can use different units which can be chosen as per user demand.

5.4.2.1 Accessing the Interface

Click **Parameter Unit** in the **Para.** selection to access the **Parameter Unit** setting interface. See Figure 5-13.




Figure 5-13 Setting Parameter Unit

5.4.2.2 Selecting Unit System

Click the **Select unit system** dropdown list and select a unit system for the parameters among the 7 unit systems (Custom, China, International, Britain, Canada, USA and Netherlands).

NOTE


- When selecting different unit standards, the corresponding unit list and unit option will be displayed differently.
- If another option is selected except the **Custom**, then the unit of each parameter can only be browsed.

5.4.2.3 Customizing Parameter Unit

1. Select **Custom** from the dropdown list of **Select unit system**.



- 2. Click the parameter, of which the unit is to be set, from the parameter list (such as **WBC**).
- 3. Select a new parameter unit from the Unit Options list.



4. Click **Apply** or **OK** to save the configuration.

NOTE

- For parameters in the same group, if the unit of any parameter changes, the units of the other parameters change accordingly. (In the list, parameters will be sorted by group; the first parameter will be displayed in black and the other parameters in the same group will be displayed in grey.)
- If the parameters units change, the display format of the list data will change accordingly.

5.4.2.4 Retrieving Defaults

When setting the **Custom** unit system, if you click **Default**, the unit of the parameters can be restored to the initial default values.

5.4.3 Ref. Range

The reference range based on various normal groups can be set for the analyzer in the actual practice. If the analysis result of a sample is beyond the reference range, it will be regarded as

clinically abnormal. The Ref. Range interface is where you view and set the high and low limits for your patients. The analyzer flags any parameter value above (\uparrow or **H**) or below (\downarrow or **L**) these limits.

This analyzer divides the patients into 4 demographic groups: **General**, **Man**, **Woman** and **Child**. You can also customize other groups. The recommended limits are for reference only. To avoid misleading parameter flags, be sure to set the patient limits according to the characteristics of local population.

5.4.3.1 Accessing the Interface

Click **Ref. Group** in the **Para.** selection to access the reference group settings interface. See Figure 5-14.




Figure 5-14 Ref. Range

5.4.3.2 Copying a Ref. Group

Select a reference group and click **Copy**, and a new reference group with everything the same except the name of the reference group will be added to the system and a screen as shown in Figure 5-15 will pop up.





Figure 5-15 Copying a Ref. Group

You can edit the new reference group. Save and close the screen, and then the copied reference group will be shown in the reference group list.



NOTE

The reference group name entered is not allowed to be empty nor the same as the existing ones.

5.4.3.3 Adding a New Ref. Group

If the built-in reference groups cannot meet the actual demand, you can add new ones and manually enter the information such as reference ranges for each parameter, names and genders. The procedures are shown as below:

1. Click **New**, and a screen for adding a new reference group will pop up. See Figure 5-16.




Figure 5-16 Adding a New Ref. Group

2. Complete the entries for each parameter with reference to the parameter description in Table 5-2.

Table 5-2 Description of Ref. Group parameters

Parameter	Meaning	Operation	
Ref. Group	Name of the new reference group.	Click the edit box and enter the information using the soft keyboard. English characters and numbers are allowed to be entered, while special characters are not.	
		NOTE	
		The reference group name entered is not allowed to be empty nor the same as the existing ones.	
Lower Limit of Age	Lower limit of age of the reference group.	Enter an integer value in the textbox and select the age unit from the drop list on the right. (Year, Month, Week, Day or Hour)	
		NOTE	
		The Lower Limit of Age must be smaller than the Upper Limit of Age.	
Upper Limit of Age	Upper limit of age of the reference group.	Enter an integer value in the edit box and select the age unit from the drop list on the right: Year , Month , Week , Day or Hour .	
		NOTE	
		The Upper Limit of Age must be greater than the Lower Limit of Age.	

Parameter	Meaning	Operation
Gender	Gender of the reference group.	Select Man, Woman, Not defined from the dropdown list. The default setting is empty.
Lower Limit (of parameter)	Lower limit of parameters of the reference group. If the test result is lower than this value, it would be regarded as clinically abnormal.	Click the Lower Limit cell which corresponds to the parameter and enter a new value. NOTE The Lower Limit must be smaller than the Upper Limit .
Upper Limit (of parameter)	Upper limit of parameters of the reference group If the test result is higher than this value, it would be regarded as clinically abnormal.	Click the Lower Limit cell which corresponds to the parameter and enter a new value. NOTE The Upper Limit must be greater than the Lower Limit .

- 3. Click **Save** to save the settings.
- 4. Click **Close** to exit the interface.

5.4.3.4 Editing a Ref. Group

You can modify the reference range of the parameters according to actual needs and set suitable reference intervals (age range, gender, etc.).

The procedures are shown as below:

1. Select the Ref. group to be set, and click **Edit** to enter the interface as shown in Figure 5-17.




Figure 5-17 Editing a Ref. Group

2. Refer to Table 5-2 for the description of the parameters to finish the editing.

NOTE

- For the built-in reference group, you can modify the upper limit and lower limit of the parameters, but not its name, the upper limit and lower limit of age as well as gender.
- Click Set as default to restore the setting of the selected reference group to the default value.
- Non-built-in reference group (which is added by user) cannot restore defaults.
- 3. Click Save to save the modification.
- 4. Click Close to exit.

5.4.3.5 Deleting a Ref. Group

Click **Delete**, and select **Yes** in the pop-up dialog box to delete the selected customized reference group.

NOTE

Built-in reference group can not be deleted.

5.4.3.6 Setting Default Ref. Group

When you pre-enter patient information in the **Sample Analysis** interface, the **Ref. Group** displayed by default is the default reference group.

The default setting is **General**. You can change it as required. Select a reference group and click **Set** as **default** to set the selected reference group as the default group.

As shown in Figure 5-18, the reference group with a check mark in its **Default** column is a default reference group.




Figure 5-18 Setting Default Ref. Group

5.4.3.7 Automatically match the customized reference group according to age and gender

If Automatically match the customized reference group according to age and gender is checked, the customized reference group will be automatically assigned patients by the system according to their age and gender when the patient information is entered. If it fails to find a matching customized reference group for a patient, the patient will be assigned to the built-in reference group.

When the system automatically matches the reference group according to age and gender, the rules listed in Table 5-3 shall be followed.

· · · · · · · · · · · · · · · · · · ·			
Automatically match the customized reference group according to age and gender	Customized reference group	Match the reference group	
Unchecked	N/A	Built-in reference group	
Checked	Does not exist	Built-in reference group	
Checked	Created	Preferentially match the customized reference group	

Table 5-3 Rules for Matching the Reference Group

NOTE

When the customized ref. groups are used to match the reference group, the matching will be performed from top down according to the customized ref. groups displayed in the screen.

5.4.4 Microscopic Exam. Settings

You can perform the microscopic exam. settings, including adding, editing and deleting as per the actual demand.

NOTE

The operations of adding, editing and deleting do not affect the sample record in which the microscopic examination results have been entered and saved. Such operations are only valid for the record in which the microscopic examination results have not been saved, and the samples analyzed after the setting operations.

5.4.4.1 Accessing the Interface

Click **Microscopic Exam.** in the **Para.** selection to access the microscopic examination setting interface. See Figure 5-19.




Figure 5-19 Microscopic Exam. Settings


5.4.4.2 Adding a New Microscopic Exam. Parameter

Do as follows to add a new microscopic exam. parameter.

1. Click **New** in the **Microscopic Exam. Settings** interface.

A dialog box will pop up as shown in Figure 5-20.

Figure 5-20 Adding a New Microscopic Exam. Parameter




- 2. Input the parameter name and its coding system in the corresponding textboxes.
 - The Parameter Name can not be empty.
 - The **Code System** is the code ID of the parameter. It is for LIS transmission only when the parameter is transmitted to the LIS. You may not input the value if it is not needed.
- 3. Click OK.

The name of the new parameter will be displayed in the microscopic exam. parameter list.

5.4.4.3 Editing a Microscopic Exam. Parameter

Select a parameter name from the list and click **Edit** to modify it. See Figure 5-21.

Figure 5-21 Editing a Microscopic Exam. Parameter



5.4.4.4 Deleting a Microscopic Exam.Parameter

Select a parameter name from the list, click the **Delete** button and then click **Yes** in the popup dialog box to delete this parameter.

Are you sure to remove the Neutrophilic segmented granulocyte?

Yes

No

Figure 5-22 Deleting a Microscopic Exam.Parameter

5.4.5 Customized Parameters

Except for this analyzer's analysis parameters, parameters collected from other testing instruments or via manual testing by the user are customized parameters. You can set customized parameters so they can be printed together with this analyzer's analysis parameter details on the **Hematology Analysis Report**.

This analyzer's default customized parameters include: Blood Type, RH Blood Group, ESR, C-reactive Protein and Reticulocyte. You can set the unit and reference range of default customized parameters as well as add and set customized parameters.

5.4.5.1 Accessing the interface

Click **Custom Para.** in the **Para.** selection to access the customized parameters setting interface. See Figure 5-1.

Custom Para. No. Parameter Name Unit 1 Blood Type RH Blood Group 2 3 **ESR** 4 C-reactive Protein 5 Reticulocyte Edit New Delete Close

Figure 5-1 Customized Parameter Settings

5.4.5.2 Adding a Customized Parameter

1. Click New,

A dialog box will pop up as shown in Figure 5-2.




Figure 5-2 Adding a Customized Parameter

- 2. Click the textboxes of **Parameter Name** and **Unit** respectively, and enter the name and unit of the customized parameter.
- 3. Click corresponding cells of the **Upper Limit** and **Lower Limit** of the reference group, and input values.

You can also customize the reference group according to the actual situation. For details, see *5.4.3 Ref.* Range.

4. Click OK.

The added parameter will be displayed in the customized parameter list.

5.4.5.3 Editing a Customized Parameter

You can set the unit and reference range of customized parameters. Detailed steps are as follows:

1. Select the customized parameter to be edited, and click **Edit**.

A dialog box will pop up as shown in Figure 5-23.




Figure 5-23 Editing a Customized Parameter

- 2. Click the textboxes of **Parameter Name** and **Unit** respectively, and modify the name and unit of the customized parameter.
- 3. Click corresponding cells of the **Upper Limit** and **Lower Limit** of the reference group, and input values.

You can also customize the reference group according to the actual situation. For details, see **5.4.3** *Ref.* Range.

4. Click Apply or OK.

5.4.5.4 Deleting a Customized Parameter

Select a customized parameter, and click on **Delete**. Then, the parameter and its corresponding reference group can be deleted.

5.5 Meterage Settings

5.5.1 Gain Settings

You can adjust each digital pot at the **Gain Settings** interface. It is not recommended to adjust gains frequently.

Click Gain Settings in the Meterage selection to access the gain setting interface. See Figure 5-24.

Gain Settings Current Value Adjustment Rate ltem **WBC** % 35 100 RBC % 135 100 **HGB Current Value** 74 HGB Blank Voltage 0.00 Apply OK Cancel

Figure 5-24 Gain Settings

NOTE

New value of the gain adjustment = Current Value × Adjustment Rate.

Setting the WBC gain

The WBC gain here is in Whole Blood Mode.

Setting method I: click the current value of the WBC and enter the new value.

Setting method II: click the **Adjustment Rate** cell of the WBC and enter the adjustment rate of the new value relative to the current value.

Setting the WBC gain

RBC channel gain.

Setting method I: click the current value of the RBC and enter the new value.

Setting method II: click the **Adjustment Rate** cell of the RBC and enter the adjustment rate of the new value relative to the current value.

Setting the HGB gain

Current digital circuit gain. The purpose for adjusting the HGB channel gain is to change the HGB background voltage.

You can enter the value directly in the **HGB Current Value** textbox or click the adjusting button to adjust the HGB gain.

Setting the HGB Blank Voltage

The background voltage derived from HGB gain cannot be modified. **HGB Background Voltage** can be adjusted within the specified range (4.2~4.8V) by modifying **HGB Current Value**.

5.5.2 Flag

When the test result meets the requirement of the Flag rules, the corresponding Flag will be displayed on the screen. You can edit the Flag rules as per the actual demand and relevant lab procedures.

Accessing the Interface

Click **Flag** in the **Meterage** selection to access the Flag rules setting interface. See Figure 5-25.




Figure 5-25 Flag

Setting Flag Rules

You can select the name of the **Flag** in the Flag interface, then click **Edit** to modify the rules in the popup dialog box. See Figure 5-26.




Figure 5-26 Setting Flag Rules

Restoring Defaults


Click **Set as default** to restore the parameter to the default value.

5.6 Communication

5.6.1 Host Network Settings

On the host communication screen, you can set the network information of the analyzer to enable its network connection.

Click **Host Communication** in the **Communicate** selection to access the host network setting interface. See Figure 5-27.



Refer to Table 5-4 for the description of relevant parameters.

Table 5-4 Description of Host Communication Setting Parameters

Parameter	Meaning	Operation
Obtain an IP address automatically	The host gets the IP address dynamically from a DHCP server or a PPP dial-up network access server.	Please choose according to the actual situation.
	This option is not applicable for the dial-up connection of SLIP server.	

Parameter	Meaning	Operation
Use the following address:	Specify the host to use the manually set IP address. If this option is selected, you need to set: IP address The IP address obtained from the network administrator or Internet service provider. Subnet mask The subnet mask obtained from the network administrator or Internet service provider. Default gateway The IP address of the default gateway; the router's IP address for connecting the independent IP network segment.	Obtain the IP address, subnet mask and default gateway of the host from the network administrator or Internet service provider.
Obtain DNS server address automatically	Automatically obtain the IP address of the Domain Name Server.	Please choose according to the actual situation.
Use the following DNS server addresses:	Specify the IP address of the DNS server of the host. Preferred DNS server The IP address of preferred or primary DNS servers. Alternate DNS server (Optional) The IP address of alternative or secondary DNS servers of the host. This server will be used if the specified IP address of the Preferred DNS server is not available or if the DNS name cannot be resolved as the IP address of the DNS server which the host has inquired.	Obtain the IP address of DNS server from the network administrator or Internet service provider.

NOTE

You can click **Details** to check the network information of the analyzer, including physical address, IP address, subnet mask, default gateway, DNS server, etc.

5.6.2 LIS Communication

In the **LIS Communication** interface, You can set the communication between the system and the LIS, including network settings, protocol settings and transmission mode.

Click **LIS Communication** in the **Communication** selection to access the Laboratory Information System (LIS) communication setting interface. See Figure 5-28.




Figure 5-28 Setting LIS Communication

Refer to Table 5-5 for the description of relevant parameters.

Table 5-5 Description of LIS Communication Setting Parameters

Parameter		Meaning	Operation
Network Settings	IP address	The IP Address of the LIS.	Please set it according to the actual situation.
	Port	The port of the LIS. The default value is 5600 .	Please set it according to the actual situation.
			An integer between 1025 and 65535 can be entered.
			NOTE
			If the analyzer is disconnected with the LIS, click the Reconnect button to connect the LIS again.
Protocol Settings	Protocol Type	Type of the protocol used for the communication between the system and the LIS. The default value is HL7 .	N/A
	Protocol Version	Version of the protocol used for the communication between the system and the LIS. The default value is 2.3.1 .	N/A

Parameter		Meaning	Operation
	Communication Acknowledgement	If checked, the communication between the system and the LIS is successful when the ACK response from the LIS is received within the duration of ACK timeout ; no response received indicates communication failure.	Please choose according to the actual situation.
		If unchecked, the communication between the system and the LIS shall be considered successful no matter the ACK response from the LIS is received or not.	
		NOTE	
		The system will send the next message continuously no matter the communication is successful or not.	
	ACK timeout	Timeout duration of the ACK response. The default value is 10 seconds, that is, the	Click or or directly input in the textbox.
		communication will be considered failed if the system	An integer between 1 and 120 can be entered.
		receives no ACK response within 10 seconds.	Unit: Second (sec.)
		within 10 seconds.	NOTE
			The parameter is valid only when the Communication Acknowledgement is checked.
Transmission Mode	Auto-communication	If checked, the system will automatically upload the result to the LIS upon the completion of the analysis.	Please choose according to the actual situation.
		If unchecked, the result of analysis will not be automatically uploaded.	
		NOTE	
		If the Bidiectional LIS/HIS Communication is checked, this parameter will be checked automatically.	

Parameter		Meaning	Operation
	Bidiectional LIS/HIS Communication	If checked, the system will automatically obtain the sample/patient information from LIS/HIS after the Sample ID is inputted or scanned, and automatically upload the result to the LIS upon the completion of the analysis.	Please choose according to the actual situation.
		If unchecked, the software system will not execute any operations.	
Histogram Transmission Method		The methods for transmitting the histogram to the LIS when the result is transmitted by the system, including:	Please choose according to the actual situation.
		Not transmit	
		Do not transmit the histogram to the LIS.	
		Bitmap	
		Transmit the histogram to the LIS in the format of screen display.	
		Transmitting bitmap for printing	
		The histogram is transmitted by the system to the LIS in the format of a printed report.	

5.7 User Management

After logging in the system, the administrator has the access to set the account information of general users and other administrators; common users can only browse the user list and change their own passwords.

5.7.1 Accessing the Interface

Click **User** in the **Setup** interface to access the user management interface as shown in Figure 5-29.




Figure 5-29 User management

5.7.2 Creating a User

Click **New** to set the account information of a new user in the popup interface, including username, first and last name, password, user group and remarks, etc. See Figure 5-30.




Figure 5-30 Creating a user

NOTE

User Group includes **Common User** and **Administrator**. Users are assigned different access levels according to the user group they belong to.

Click **OK** after the setting is complete. The information of the new user will be shown in the user list.

5.7.3 Editing a User

Select the user to be edited and click **Edit** to modify the name and user group.




Figure 5-31 Editing a User

5.7.4 Deleting a User

Select the user to be deleted and click **Delete**, and then select **OK** in the pop-up dialog box to delete the user.




The administrator cannot delete his/her own information.

5.7.5 Setting the Default User


Select a user and click **Set as default user** to set this user as the default user.

After the cleaning is completed, the following message box will pop up.



After it is set successfully, the default user name will be displayed in the login box next time and you only needs to enter the corresponding password. See Figure 5-32.


Figure 5-32 Login after Setting the Default User



5.7.6 Changing Password

Click **Change Password**, enter the old password and new password of the user and confirm the new password in the popup dialog box, then click **OK**.

Figure 5-33 Changing Password



NOTE

You can only change his/her own password and cannot change the password of other users.

5.7.7 Resetting Password

If the user forgets the password or the password is required to be reset due to other reasons, please click **Reset Password** to reset the password of the selected user to the initial password. The reset password is the same as the user name.

Figure 5-34 shows that the password is successfully reset.

Figure 5-34 Resetting Password



NOTE

The administrator is allowed to reset the password of all administrators and general users; general users do not have the access to reset the password.

5.8 Print Settings

Click **Print Settings** in the **Setup** interface for relevant print settings, including the default printer, template, report, copies and margins, etc.




Figure 5-35 Print Settings

Printer Settings

You can set the printer and driver of the system in the **Printer** combo box. See Figure 5-36.

Figure 5-36 Printer Settings



Printer Driver

The analyzer uses the built-in thermal printer by default. If the **Check automatically** is selected, you can select a different printer in the **Printer** list.

Printer

The analyzer uses the built-in thermal printer by default. Select a printer to be used from the dropdown list as required. If the dropdown list is blank, it indicates that no printer has been installed for the operating system. In this case, install a printer, and then perform the relevant settings and printing operations.

Report Settings

You can set relevant parameters of the report in the Report Settings combo box. See Figure 5-37.

Report Settings

Report Title lematology Analysis Report

Copies 1

Report Title

Enter the title of the report in the **Report Title** textbox. The default setting is **Hematology Analysis Report**.

Copies

You can enter the number of copies to be printed for a report in the **Copies** textbox according to the actual demand. Click to increase the number of copies and click to decrease the number of copies or enter the number of copies in the edit box directly. Range of the copies is between 1 and 100 and the default value is **1**.

Format Settings

Report type and template of prints can be set in the Format Settings combo box. See Figure 5-38.

Format Settings

Report Type Report

Template A4-Portrait-Paramete

Paper Type A4

Paper size 210*297 mm

Refresh Import Set as default

Figure 5-38 Format Settings

Selecting Report Type

Select the format type to be set from the dropdown list of the **Report Type**. The default setting is **Report**.

Selecting Template

Select the template to be set from the dropdown list of the **Template**. The default value is **A4-Portrait-Parameters-Microscopic Exam.**.

After the selection is completed, paper type and size of the template will be shown at the bottom of the list.

Refresh

Click **Refresh** to refresh the format list after the customization by the administrator.

Importing template

Click **Import** to import the report template from the USB flash disk.

Setting the Default Template

Select the report template according to the actual demand and click **Set as Default** to set the current template as the default template.

Auto Settings

Autoprint

The default setting is **Off**, which means the report should be printed manually after the results are obtained.

If it is set to **On**, the system will automatically print the report of the sample as per the current report template once the counting results are obtained.

NOTE

- If **Print after validation** is checked, the autoprint function becomes invalid.
- Auto print is not applicable for the background results.
- Auto print after validation

It's unchecked by default, which means the system can print the report automatically without validation.

If it's checked, the report will be printed automatically after it's been validated instead of being printed right after the results are obtained each time.

NOTE

The parameter is valid only when the **Autoprint** is set to **On**.

Auto validate when printing

It's unchecked by default, which means the report will not be automatically validated by the system at the time of printing.

If it's checked, the report will be automatically validated and printed by the system at the time of printing.

Print after validation

It's unchecked by default, which means the report can be printed without validation.

If it's checked, the report can be printed only after validation and autoprint is unexecutable.

Printing Options

Print Flag

It's checked by default, which means the flag information will be printed in the report. If it's not checked, it will not be printed.

Print Ref. Range

It's checked by default, which means the reference range of the parameter will be shown in the printed report; If it's unchecked, the results alone, rather than reference range, will be shown in the printed report and the reference range will not.

Print Suspicious Flag

It's unchecked by default, which means the suspicious flag "?" will not be shown in the printed report; if it's checked, such flag can be shown.

Print Ref. Range Flags

It's checked by default, which means the printed report can show the ref. range flag (such as ↑ or ↓); if it's unchecked, such a flag will not be shown.

Print Microscopic Exam. Para.

It's checked by default, which means the result of Microscopic Exam. Parameters will be printed in the report. If it's not checked, it will not be printed.

Print result edited flags

It's unchecked by default, which means the mark for the edited results will not be shown in the printed report.

If checked, the mark (**M** or **m**) for the edited results will be shown in the printed report if the parameters have been modified.

Two reports in one page (half of A4)

It's unchecked by default. If this is checked, the default template size in **Format Settings** is half an A4 page (e.g., **A4_Half-Portrait-Parameters**), so two reports can be printed in one piece of A4 paper.

NOTE

When Autoprint is On, a page remains to be printed with one report.

Update blank test time before be printed

It's unchecked by default, which means the blank test time will not be processed by the system. If it's checked, the **Delivery Time** will be automatically updated as the **Run Time** by the system at the time of printing.

Print as black and white

It's unchecked by default, whick means the report will be printed according to the default settings of the printer.

If it's checked, the resport will be printed as black and white.

5.9 Auxiliary Settings

Click **Auxiliary Settings** in the **Setup** interface to access the **Auxiliary Settings** interface. See Figure 5-39.




Figure 5-39 Auxiliary Settings

The administrator is allowed to set the following functions in the Auxiliary Settings interface:

- Sample numbering rules
- Startup sample IP and mode
- Predilute
- Other

Sample Numbering Rules

Set the sample ID entry rules.

Sample ID Entry Method

Click the dropdown list of the **Sample ID Entry Method** and select the entry method of the sample ID from the following options.

- Auto increment (default setting)
- Manual entry

Prefix Length

When **Auto Increment** is selected as the **Sample ID entry method**, you can add a prefix to a certain batch of samples for identification.

Enter the prefix length ranging from 0 to 24 (e.g. 2) of the sample ID in the **Prefix Length** textbox. The prefix length will applied to all sample IDs after the setting is saved.

Startup sample IP and mode

Set the sample ID and measurement mode for the next sample after startup.

Next Sample ID and mode after startup

The sample ID set by the user will be used after the next startup when the specified sample ID is entered into the textbox.

NOTE

If the **Effective tomorrow** is checked, the modification of the next sample ID and mode after startup will become effective on the next day.

Continue using the sample ID and mode before the last shutdown
 If checked, the system will by default add 1 to the last sample ID analyzed before shutdown as the next sample ID after startup.

Predilute

Set if you wish to see a popup dialog box when you perform the Predilute counting.

- Ask for confirmation (default setting): in the **Predilute** mode, when you press the aspirate key to start the analysis, a dialog box will pop up to remind you that the ongoing analysis is for Predilute counting.
- Do not ask for confirmation: the dialog box for confirming the Predilute counting will not pop up.

Other

Show Result Edited Flags

It's unchecked by default, which means the edited results are marked with an $\bf M$ at the end, while the corresponding results with manual modifications are marked with an $\bf m$ at the end. $\bf M$ or $\bf m$ is displayed between the result data and the parameter unit by default.

If unchecked, the edited result will not be marked with an **M** or **m**.

Suspicious Flag

A single character (an English letter only) can be re-entered in the textbox as a suspicious flag. The default value is **?.**

● Ref. Range Flags

ਮੁੰਦੂ ਅੰਗਰ ਜ਼ਿਲ੍ਹਾਂ ਜ਼ਿਲ੍ਹ

6 Daily Operations

6.1 Introduction

This chapter introduces the daily operations from the startup to the shutdown of the analyzer.

A flow chart indicating the common daily operation process is presented below.




Figure 6-1 Daily Operations Procedure

6.2 Pre-operation Preparation



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.



WARNING

- Be sure to dispose of reagents, waste, samples, consumables, etc. according to local legislations and regulations.
- The reagents are irritating to eyes, skin and mucosa. Wear proper personal protective
 equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when
 handling them in the laboratory.
- If the reagents accidentally spill on the skin, wash them off with plenty of water and if necessary, go see a doctor; if the reagents accidentally spill into the eyes, wash them off with plenty of water and immediately go see a doctor.
- Keep your clothes, hairs and hands away from the moving parts to avoid injury.
- The sample probe tip is sharp and may contain biohazardous materials. Exercise caution to avoid contact with the probe when working around it.

NOTE

- You should only use the MR-specified reagents. Store and use thereagents as specified in instructions for use of the reagents.
- Check if the reagents are connected correctly before using the analyzer.
- After long-distance transportation, the reagent must be allowed to settle for more than one day before use.
- Be sure to use clean K₂EDTA vacutainer blood collection tubes with anticoagulant, fused silica glass/plastic test tubes, centrifugal tubes and borosilicate glass capillary tubes.
- Be sure to use the MR-specified disposable products including vacutainer blood collection tube, vacutainer blood collection tubes with anticoagulant and capillary tubes etc.

Perform the following checks before turning on the analyzer.

Waste container

Check and make sure the waste container is empty.

- Fluidic tubing and power connections
 - > Check and make sure the reagents and waste tubing are properly connected and not bent.
 - Check and make sure the power cord of the analyzer is properly plugged into the power outlet.
- Printer (Optional)

Check and make sure enough paper is installed.

Network Cable (Optional)

Check and make sure the network cable is properly connected to the analyzer.

6.3 Startup

This section introduces the operations related to the startup of the analyzer.

NOTE

- If you failed to start the analyzer continuously, please contact MR customer service department or your local agent immediately.
- After startup, please make sure the data/time displayed on the screen is correct.

1. Place the power switch at the back of the analyzer in the [I] position.


The power indicator light will be on.

2. Check the indicator light on the analyzer.

If the indicator light is on, it indicates the analyzer has been started up. The analyzer will perform self-test and initialization in sequence. The whole process will last for 4 to 10 minutes. (Time needed for initializing the fluidic systems depends on how the analyzer was previously shut down.)

3. Enter the correct user name and password in the Login message box. See Figure 6-2.

Figure 6-2 Login



The initial user name and password of administrator are **admin**, which was set by service engineer.

1 to 12 digits of numeric characters can be entered for the user name and the password. No Chinese character is allowed.

4. Click to enter the user interface.

The system will display the **Sample Analysis** screen by default and display the test result of the background when the analyzer is started.

NOTE

- The Sample Analysis interface will be displayed by default. And the background checked will
 be displayed on the screen. The background test is designed for detecting particle interference
 and electrical interference.
- For the background Ref. Range of each parameter, please see A.4.2 Normal Background.
- The sample ID for the background test is background.
- If the background results exceed the Ref. Range for the first time during fluidics initialization, then the analyzer will run the background test one more time.
- Running a test when there is a **Background abnormal**, you would obtain an unreliable testing result.
- If any error is detected during initialization (e.g. the background results exceed the Ref. Range), the analyzer will activate the alarm. For details, see *13 Troubleshooting*.
- To lock or switch a user, click on the menu screen and click **Yes** on the pop-up dialog box. The system will return to the login dialog box. Enter the user name and password, click then you can log in again or log in the software interface with another user identity.

6.4 Daily Quality Control

To ensure reliable analysis results, conduct daily QC analysis on the analyzer before running samples. For details, see **9 Quality Control**.

6.5 Sample Collection and Handling



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.



WARNING

Do not touch the patients' blood sample directly.



CAUTION

- Do not re-use such disposable products as collection tubes, test tubes, capillary tubes, etc.
- Prepare the samples as per the procedures recommended by the reagent manufacturer.

NOTE

- Be sure to use clean K₂EDTA vacutainer blood collection tubes with anticoagulant, fused silica glass/plastic test tubes, centrifugal tubes and borosilicate glass capillary tubes.
- Be sure to use the MR-specified disposable products including vacutainer blood collection tube, vacutainer blood collection tubes with anticoagulant and capillary tubes etc.
- For the whole blood samples to be used for WBC classification or PLT count, store them at room temperature and run them within 8 hours after collection.
- If you do not need the PLT, MCV and WBC differential results, you can store the samples in a refrigerator (2°C - 8°C) for 24 hours. You need to warm the keep samples at room temperature for at least 30 minutes before running them.
- Be sure to shake any sample that has been prepared for a while before running it.

6.5.1 Running the Venous Whole Blood Samples

The procedure for preparing capillary whole blood sample is as follows:

- 1. Use clean K₂EDTA (1.5~2.2mg/mL) vacutainer blood collection tubes with anticoagulant to collect venous blood samples.
- 2. Mix the venous blood with the anticoagulant well in the tube immediately.



CAUTION

For vacutainer blood collection tube (Φ 12X75, cap excluded), please make sure the volume of the whole blood sample is not less than 0.5mL.

6.5.2 Capillary Whole Blood Samples

Collect the capillary whole blood sample with a vacuum blood collection tube specified by the manufacturer.



CAUTION

To ensure the accuracy of the analysis, make sure the volume of the capillary whole blood sample is not less than $100\mu L$.

NOTE

Run the capillary whole blood sample within 3 minutes to 2 hours after its collection.

6.5.3 Prediluted Samples

The procedure for preparing prediluted sample is as follows:


1. Click the on the top left corner and enter the menu screen as shown in Figure 6-3




Figure 6-3 Menu Screen

2. Click the Add Diluent icon.

A prompt box will pop up on the screen as shown below.



3. Take a clean centrifugal tube, uncap it and present it to the sample probe in a manner as shown in the following picture in which the probe tip is vertically in contact with the bottom of the tube so as to avoid bubbles, liquid attached to the inner wall or spatter.



- Press the aspirate key and add the diluent (180μL at a time).
 After the diluent is added and you hear a beep, you can remove the centrifugal tube.
- 5. If more portions of diluent are needed, repeat steps 3~4.
- 6. Add 20µL of blood to the diluent, close the tube cap and shake the tube to mix the sample.
- 7. After the prediluted sample is prepared, click **Cancel** to exit dispensing the diluent.

NOTE

- You can also dispense 180μL of diluent by pipette into the tube.
- Be sure to keep dust from the prepared diluent.
- Be sure to run the prediluted samples within 30 minutes after the mixing.
- Be sure to mix any sample that has been prepared for a while before running it.
- Be sure to evaluate predilute stability based on your laboratory's sample population and sample collection techniques or methods.

6.6 Sample Analysis

After the sample is prepared, you can perform the operations for sample analysis.

For details, see 6.6 Sample Analysis.

6.7 Shutdown



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.



WARNING

The sample probe is sharp and potentially biohazardous. Do not turn on the analyzer immediately after its shutdown.



CAUTION

Wait at least 10 seconds before power-on to avoid damage to the machine.


NOTE

- To ensure stable analyzer performance and accurate analysis results, be sure to perform the Shutdown procedure to shut down the analyzer after it has been running continuously for 24 hours.
- When the analyzer is running or performing other fluidics sequence, do not force shutdown the analyzer.
- If any error is detected during shutdown procedure, the analyzer will return to the status before
 the shutdown procedure is performed, and then activate the alarm. See 13 Troubleshooting for
 details of removing the error.
- Be sure to shut down the analyzer in strict accordance with the instruction below.

Procedures for shutting down the analyzer are as follows:

1. Click the button on the menu screen.

A dialog box will pop up as shown below.



2. Click Yes.

The system starts to execute the shutdown sequence and a message box pops up showing the procedures for cleanser maintenance.

3. Follow the instructions and set the cleanser under the sample probe, and press the aspirate key on the analyzer or click **Aspirate** to run the cleanser aspiration.

Upon the completion of cleanser maintenance, a message displayed on the screen indicates the cleanser maintenance is completed.

Shutdown done. Please power off the analyzer!

- 4. Place the [O/I] switch at the back of the main unit in the [O] position.
- 5. After shutdown, empty the waste in the waste container, and dispose of it.



WARNING

Be sure to dispose of reagents, waste, samples, consumables, etc. according to local legislations and regulations.

7 Sample Analysis

7.1 Introduction

Sample analysis is the most important function of the auto hematology analyzer. You can get the blood cell count, HGB concentration and the 3-part classification counting results of the white blood cells by performing the sample analysis.

The summary of sample analysis procedures is as follows:

- 1. Entering the sample information.
- 2. Running the samples.
- 3. Processing the analysis results.

7.2 Interface Introduction

The **Sample Analysis** interface is the main interface of the analyzer (Figure 7-1). You can complete the operations such as entering the sample information, performing sample analysis, reviewing/printing analysis results in the **Sample Analysis** interface.



Figure 7-1 Sample analysis interface

Related descriptions:

Function buttons

You can perform operations such as setting the mode for the samples, pre-entering information, reviewing previous/next records and printing. Click and view all function buttons. See section **7.6** *Functions of the Buttons*.

Patient information area

It displays the patient information corresponding to the current sample.

Analysis results area

It displays the analysis results of the sample, including the parameter results, Flags and histograms. They system displays the analysis results of the most recent run by default.

Parameter Results

This list displays the analysis results of all the parameters of the samples.

You can compare the values in the Result column with the corresponding Ref. Range. If the values are within the reference range, it means that they are normal. If not, it indicates that the sample may be abnormal and the corresponding symbols will be displayed in the Flag column.

- WBC Message Displays the alert message regarding the WBC.
- > RBC Message

Displays the alert message regarding the RBC.

PLT Message

Displays the alert message regarding the platelet.

➤ WBC

RBC distribution histogram

➤ RBC

RBC distribution histogram

▶ PLT

Platelet distribution histogram.

Information of the next sample

It displays the sample ID and analysis mode of the next sample.

7.3 Entering Sample Information

You can enter the worklist information of the samples to be tested before the analysis.

NOTE

You can also enter sample/patient information after the sample analysis is completed. For details, please refer to *8 Result Review*.

Specific steps are shown below:

1. Click the Pre-entry button in the function button area.

The interface as shown in Figure 7-2 will pop up on the screen.




Figure 7-2 Pre-entering Patient Information

2. Enter patient information with reference to the parameter description in Table 7-1.

Table 7-1 Parameter Description

Parameter	Meaning	Operation
First Name	First name of patient.	Enter in the textbox directly.
Last Name	Last name of patient.	Enter in the textbox directly.
Patient Type	Type of patient. Values: Inpatient Physical Examination STAT Outpatient	Select from the dropdown list.

Parameter	Meaning	Operation
Sample Type	Type of sample for microscopic examination. Values: • Venous blood • Capillary • Cord blood • Blood	Click the Sample Type dropdown list box and select the type of sample for microscopic examination.
Med Rec. No.	Med Rec. No. of patient.	Enter in the textbox directly.
Gender	Gender of patient. Values: • Male • Female • Not defined	Select from the dropdown list.
Birthday	Birthday of a patient.	 Select from the date control. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd, you should input the data in the sepuence of year, month, and date. Click or select the date or click the textbox to input them directly. Click to clear the data and input again.
Age	Age of a patient.	Select the unit of age from the dropdown list (Year, Month, Day or Hour) and enter the age of the patient in the textbox before the age unit. NOTE If the Birthday is set, the age will be displayed automatically.
Ref. Group	Reference group of the sample under analysis. The result is judged according to the reference range of the reference group and the result beyond the normal range will be flagged.	Select from the dropdown list. NOTE • If the Automatically match the customized reference group according to age and gender is set, gender and age of a patient will automatically match the reference group according to the corresponding relationship (No matter the reference group is selected or not). • Refer to 5.4.3 Ref. Range for the setting of the reference group and range.
Department	Department receiving the patient.	Select from the dropdown list.
Area	Ward area of patient.	Select from the dropdown list or input directly.

Parameter	Meaning	Operation
Bed No.	Bed No. of inpatient.	Select from the dropdown list or input directly. NOTE The bed No. is required to be filled only for inpatients.
Sampling Time	Date and time when the sample is collected.	Click the date control for the settings. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd HH:mm, you should input the data in the sepuence of year, month, date, hour, and minute. Click or to select the date or click the textbox to enter them directly. Click to delete the current data and re-enter information. NOTE The system automatically displays the current time as sampling time. The sampling time can be no later than the current system time.
Submitter	Personnel submitting the sample.	Select from the dropdown list or input directly.
Delivery Time	Date and time when the sample is delivered.	 Click the date control for the settings. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd HH:mm, you should input the data in the sepuence of year, month, date, hour, and minute. Click or to select the date or click the textbox to enter them directly. Click to delete the current data and re-enter information. NOTE The system automatically displays the current time as sample delivery time. The delivery time can be no later than the current system time.
Remarks	Clarifications or notes.	Input in the textbox directly.

3. Click **Apply** or **OK** to save the configuration.

7.4 Running Samples



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.



WARNING

The sample probe tip is sharp and may contain biohazardous materials. Exercise caution to avoid contact with the probe when working around it.



CAUTION

- Do not re-use such disposable products as collection tubes, test tubes, capillary tubes, etc.
- Make sure that the entered sample ID and mode exactly match those of the samples to be run.

NOTE

- During aspiration, the tip of the probe should be kept at a certain distance from the bottom of the sample container, otherwise the accuracy of aspiration volume will be affected.
- Keep the tip of the probe from contacting with the wall of the test tube to avoid blood splashing.
- Proper reference range shall be selected on the **Setup** interface before analysis. Otherwise, the
 results may be flagged erroneously.
- The default system setting for Mode & ID is Venous Whole Blood (VWB).
- When the analyzer is running the samples, you can switch to **Review** interface to perform operations including browsing and exporting, etc., and you can also switch to other interfaces.
- When the analyzer is running the samples, all the functions related to the fluidics sequence are not available.

Take the following steps to perform sample analysis.

- 1. Prepare samples as instructed by 6.5 Sample Collection and Handling.
 - For details about the preparation of venous whole blood samples, see **6.5.1** Running the Venous Whole Blood Samples.
 - For details about the preparation of capillary whole blood samples, see **6.5.2** Capillary Whole Blood Samples.
 - For details about the preparation of prediluted samples, see 6.5.3 Prediluted Samples.
- 2. Shake the capped tube of sample for a homogeneous specimen.
- When the green indicator light is steady-on, click Mode & ID in Sample Analysis interface.
 A dialog box will pop up as shown below.




Figure 7-3 Mode & ID Settings

- 4. Select the blood sample mode **Venous Whole Blood** (VWB), **Capillary Whole Blood** (CWB) or **Predilute** (PD) of the sample, and input the **Sample ID**.
 - Letters, numbers and all characters that can be entered through the keyboard (including special characters) are allowed for the Sample ID. Chinese and other languages (such as Japanese, Korean, etc) are not supported.
 - The length of the entries ranges from 1 to 25 and the entries shall not be empty.
 - ➤ The last character of a sample ID must be numeric, but a string of "0" only is not an acceptable sample ID.
- 5. Click OK.
- 6. Remove the tube cap carefully and place the sample under the probe so that the probe can aspirate the well-mixed sample.
- 7. Press the aspirate key on the analyzer to start running the sample.

 The sample will be automatically aspirated by the sample probe.
- 8. When you hear a beep, remove the sample tube.
 - The analyzer will automatically run the sample and the analysis status icon and analyzer indicator is flickering in green. When the analysis is complete, the analyzer indicator returns to constantly-on green.
- 9. Repeat steps 1~8 to run the remaining samples.

7.5 Dealing with the Analysis Results

7.5.1 Automatic saving of analysis results

This analyzer automatically saves sample results. When the maximum number has been reached, the newest result will overwrite the oldest (already backed up). The maximum number of the automatically saved results is 50,000.

7.5.2 Parameter Flags

- If parameter is followed by a "↑" (H) or "↓"(L), it means the analysis result has exceeded the
 upper or lower limit of the reference range but still within the display range.
- If the parameter is followed by a "?", it means the analysis result is suspicious.
- If you see "***" instead of a result, it means the result is either invalid or beyond the display range.

NOTE

For the background test, the flags for parameters or abnormal blood cell differential and morphology are not available.

7.5.3 Flags of Abnormal Blood Cell Differential or Morphology

The analyzer will flag abnormal or suspicious WBC, RBC and PLT according to the scattergrams and histograms. The flag information is defined in the table below.

Table 7-2 Flags of abnormal blood cell differential or morphology

Flag Type		Flag information
	Abnormal	Leucocytosis
		Leucopenia
		Granulocytosis
		Granulopenia
		Lymphocytosis
WBC		Lymphopenia
		MID Increase
	Suspicious	Background/Aspiration Abn.
		WBC abnormal?
		Abnormal WBC histogram
		Abnormal WBC Channel
	Abnormal	Erythrocytosis
		Anisocytosis
		Macrocytosis
RBC/HGB		Microcytosis
		Anemia
		Hypochromia
	Suspicious	Abnor. RBC Distr.
		Dimorphologic

Flag Type		Flag information
		Iron Deficiency?
		HGB Abn/Interfere?
		RBC Clump?
		Abnormal RBC Channel
		Abnormal HGB Channel
PLT	Abnormal	Thrombocytosis
		Thrombopenia
	Suspicious	Abnor. PLT Distr.
		PLT Clump?

7.6 Functions of the Buttons

7.6.1 Previous/Next

Click **Previous**, and the screen will display the sample analysis results prior to the current one. Click **Next**, and the screen will display the sample analysis results after the current one.

7.6.2 Mode & ID

Click this button to set the sample mode and measurement mode during the sample analysis. See section **7.4** *Running Samples*.

7.6.3 Pre-entry

Click this button, and you can pre-enter the information of the sample to be tested before performing the sample analysis. See section **7.3** *Entering* Sample Information.

7.6.4 Validate/Cancel Validation

After running sample, you can click **Validate** to validate the sample. After validating, the button will replaced by **Cancel Validation**. After validating, you can not edit the sample/patient information and the result.

If the current sample has been validated, the sample validation can be canceled by clicking **Cancel Validation**. After canceling the validation, you can edit the sample/patient information and the result.

7.6.5 Print

You can click Print to print the report of the sample result.

7.6.6 Customized Parameters

You can browse and edit the customized parameters results of the selected sample in the **Sample Analysis** interface. The operation procedures are as shown below:

1. Click **Custom Para.** to enter the customized parameters setting interface as shown in Figure 7-4.

Para. Flag Value Unit Range
Blood Type
RH Blood Group
ESR
C-reactive Protein
Reticulocyte

Apply
OK Close

Figure 7-4 Customized Parameters

2. Click the cell corresponding to its Value column of the parameter, and enter the value.

If the unit and reference range of parameters have been set in the **Setup > Parameter > Custom Para.** interface, the corresponding unit and range (lower limit~upper limit) will be displayed in this tab. When both the value and range of parameters are numbers, and the number is out of the reference range, the relevant mark \uparrow or \downarrow will be displayed in the Flag column.

Please refer to **5.4.5 Customized Parameter** for customized parameters settings.

7.6.7 Microscopic Exam. Parameters

You can perform the microscopic examination settings as per the following steps.

1. Click Microscopic Exam. Para.

The microscopic examination parameters interface as shown in Figure 7-5 will pop up on the screen.




Figure 7-5 Adding a New Microscopic Exam. Parameter

2. Set the microscopic examination parameters by referring to Table 7-3.

Table 7-3 Microscopic Exam. Parameters

Parameter	Meaning	Operation
Sample Type	Type of sample for microscopic examination. • Venous blood • Capillary • Cord blood • Blood	Click the Sample Type dropdown list box and select the type of sample for microscopic examination.
Blood Type	Blood type of a patient.	Select the blood type of the patient in the Blood Type/ESR column. Click the first combo box next to the blood type, you can select from Blank , A , B , O and AB ; click the second combo box, you can select from Blank , RH+ and RH- .

Parameter	Meaning	Operation
Exam. Time	Time of microscopic examination.	 Click the date control for the settings. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd HH:mm, you should input the data in the sepuence of year, month, date, hour, and minute. Click or to select the date or click the textbox to enter them directly. Click to delete the current data and re-enter information. NOTE The Microscopic exam. time can be no later than the current system time.
Microscopic Description	Description of cells morphology.	Enter the morphology information for cells into the multi-line textbox.

7.6.8 Communication

You can transmit the current sample data (except the background sample) to the LIS/HIS system in the **Sample Analysis** interface. The operation procedures are as shown below:

- 1. Click to unfold all function buttons.
- 2. Click Comm..

7.6.9 Edit Result

NOTE

- You can not edit the results of validated samples.
- You can not edit the results of the background.

You can edit the parameter result of the selected sample as per the following steps.

- 1. Click to unfold all function buttons.
- 2. Click Edit Result.

The Edit Result dialog box will pop up on the screen as shown in Figure 7-6.

Edit Result WBC 10⁹/L 10⁹/L 4.05 **RBC** 4.46 10¹2/L PLT 322 Lym% 35.9 % HGB 101 g/L MPV 8.3 fL GRAN% 53.5 % **HCT** 38.6 % **PDW** 16.1 MID% % RDW-CV P-LCR % 10.6 14.1 30.2 RDW-SD 52.6 fL Apply OK Cancel

Figure 7-6 Editing Parameter Result

- 3. Modify the counting results of the corresponding sample parameters.
- 4. Click **Apply** or **OK** to save the changes.

If the sum of the percentage of the diff parameters is not equal to 100.00% or the WBC value is invalid after modification, the system will prompt in a message box that the entered value is invalid. Please re-enter after confirmation.

If the result of one parameter is modified, then the result of other related parameter(s) will be changed accordingly and the high or low/suspicious flags will also be updated.

NOTE


The result of the parameter that you modified manually will be flagged with an \mathbf{M} . If any parameter result is then changed due to the one that you modified manually, it will be flagged with an \mathbf{m} .

7.6.10 Delete

NOTE

- Validated samples are not allowed to be deleted.
- The common user has no access to delete the sample records.
- 1. Click to unfold all function buttons.
- 2. Click **Delete**, and then click **Yes** in the pop-up dialog box to delete the sample.

Figure 7-7 Delete Sample Records



8 Result Review


8.1 Introduction

Upon the completion of each sample analysis, the analyzer will automatically save the sample information, result data, flag messages, histograms and scattergrams to the Review Database. The sample pool of the analyzer can save up to 50,000 sample records.

In the **Review** Interface, you can browse the saved sample information, result data, flag messages, histograms and scattergrams, and can search, compare or export the saved sample information.

8.2 Interface Introduction

You can browse, search, compare, print, and export the existing results in the Review interface. Click **Review** to enter the sample review interface. See Figure 8-1.



Interface Description:

- Result list: you can see browse detailed sample records.
- Function buttons: you can perform the operations such as comparing or searching the sample results, deleting and viewing the Run Charts, exporting and printing reports.
- Direction button: If you click different direction buttons, the list will move toward the corresponding directions.
 - From left to right, it indicates in sequence: the first column, moving to the left page, moving to the right page, and the last column.
 - From top to bottom, it indicates in sequence: the first page, the previous page, the next page, and the last page.

8.3 Sample List

The review interface shows a list of the analyzed samples, which contains the sample ID, mode, status and results of various parameters and other information.

Click a sample or multiple samples in the list area, then you perform operations such as exporting in batch for the selected samples. To cancel the selection, click the selected samples again.

8.4 Functions of the Buttons

8.4.1 Validate

NOTE

After validating, you can not edit the sample/patient information and the result.

After running samples, you can validate the samples as per the following steps.

1. Click Validate.


A dialog box will pop up as shown below.



- 2. Select the sample which needs to be validated.
 - Selected Records: The sample results shown on the highlighted page.
 - > Samples on current page: Results of all the samples shown on the current page.
- 3. Click OK.

The system will prompt the validation results as shown in Figure 8-2.

Figure 8-2 Validation Results



4. Click **OK** to close the message box.

8.4.2 Cancel Validation


NOTE

After canceling the validation, you can edit the sample/patient information and the result.

You can cancel the validation of validated samples. Specific steps are shown below:

1. Click Cancel Validation.


A dialog box will pop up as shown below.



- 2. Select one or more samples to be validated.
 - > Select **Selected Records**, and the system will cancel the validation for the current highlighted sample results.
 - Select Samples on current page, and the system will cancel the validation for all the samples on the current page.
- 3. Click OK.

The system will prompt the operation results as shown in Figure 8-3.

Figure 8-3 Validation Results



4. Click **OK** to close the message box.

8.4.3 Print

Click **Print** to print the result report of the selected sample.

8.4.4 Delete

NOTE

- Validated samples are not allowed to be deleted.
- The common user has no access to delete the sample records.
- 1. Select one or several sample records to be deleted.
- 2. Click Delete.

A prompt box will pop up on the screen as shown below.

Figure 8-4 Delete Sample Records



- 3. Select one or several sample records to be deleted according to the actual situation.
 - > Selected Records: The sample results shown on the highlighted page.
 - Samples on current page: Results of all the samples shown on the current page.
- 4. Click **OK** to delete the selected record(s).

8.4.5 Export

The operator can export the sample data to the USB flash disk for backup. There are two ways of exporting the sample data: exporting selected records and exporting records of specified dates.

- Export Selected Records
 - a. Insert a USB flash disk in the USB interface on the analyzer.
 - Select records to be backed up, and click Export.
 As shown in the following figure, the export range of the system is Selected Records by default.





Figure 8-5 Export Selected Records


- Select the content to be exported according to the actual demand.
 Content available for export includes: Patient Info., Sample Info., Graphs and Flags, and Custom Para..
- d. Click OK.
- e. Select the data export path in the popup dialog box, enter the backup file name, and click **Save**.

The file will be exported to the root directory of the USB flash disk (/udisk/sda1) and named in the format of SampleInfo_yyyyMMdd_hhmmss.csv. Among which, yyyyMMdd_hhmmss means data export year, month, date, hour, minute, and second.



f. Click Save.

The system pops up a dialog box as shown below to indicate that the data export is successful.



- Export Records of the Specified Dates
 - a. Insert a USB flash disk into the USB interface on the analyzer.
 - b. Click Export.
 - c. Select **Records of the Specified Dates** and set the run date range of sample in the two date textboxes. See Figure 8-6.




Figure 8-6 Export Records of the Specified Dates

- d. Select the content to be exported according to the actual demand.
 Content available for export includes: Patient Info., Sample Info., Graphs and Flags, and Custom Para..
- g. Click OK.
- h. Select the data export path in the popup dialog box, enter the backup file name, and click Save.

The file will be exported to the root directory of the USB flash disk (/udisk/sda1) and named in the format of SampleInfo_yyyyMMdd_hhmmss.csv. Among which, yyyyMMdd_hhmmss means data export year, month, date, hour, minute, and second.

i. Click Export.

The system pops up a dialog box as shown below to indicate that the data export is successful.



8.4.6 Edit Result

NOTE

- You can not edit the results of validated samples.
- You can not edit the results of the background.

You can edit the parameter result of the selected sample as per the following steps.

1. Select a row of record from the result list and click the Edit Result button.

The Edit Result dialog box will pop up on the screen as shown in Figure 8-7.

Edit Result **WBC** 10³/uL 10⁴3/uL **RBC** 10^6/uL PLT 6.57 4.62 284 Lym% g/dL HGB MPV fL 37.4 15.1 9.0 GRAN% HCT **PDW** 59.1 46.7 15.5 MID% % RDW-CV P-LCR % 3.5 13.3 37.3 RDW-SD 57.5 fL Apply OK Cancel

Figure 8-7 Editing Parameter Result

- 2. Modify the counting results of the corresponding sample parameters.
- 3. Click **Apply** or **OK** to save the changes.

If the sum of the percentage of the diff parameters is not equal to 100.00% or the WBC value is invalid after modification, the system will prompt in a message box that the entered value is invalid. Please re-enter after confirmation.

If the result of one parameter is modified, then the result of other related parameter(s) will be changed accordingly and the high or low/suspicious flags will also be updated.


8.4.7 Patient Info.

You can browse and edit sample/patient information after the sample analysis is completed. Specific steps are shown below:

1. Click Patient Info..

The patient information page will pop up as shown in Figure 8-8.

Figure 8-8 Patient Info.



2. Enter patient information with reference to the parameter description in Table 8-1.

Table 8-1 Parameter Description

Parameter	Meaning	Operation			
Sample ID	Number of the selected sample.	It will be displayed automatically, and you can modify it manually.			
First Name	First name of patient.	. Input in the textbox directly.			
Last Name	Last name of patient.	Input in the textbox directly.			
Patient Type	Type of patient. Values: Inpatient Physical Examination STAT Outpatient	Select from the dropdown list.			

Parameter	Meaning	Operation
Sample Type	Type of selected sample. • Venous blood • Capillary • Cord blood • Blood	Select from the dropdown list.
Med Rec. No.	Medical record number of patient.	Input in the textbox directly.
Gender	Gender of patient. Values: • Male • Female • Undefined	Select from the dropdown list.
Birthday	Birthday of a patient.	 Select from the date control. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd, you should input the data in the sepuence of year, month, and date. Click or select the date or click the textbox to input them directly. Click to clear the data and input again.
Age	Age of a patient.	Select the unit of age from the dropdown list (Year , Month , Day or Hour) and enter the age of the patient in the textbox before the age unit.
Ref. Group	Reference group of the sample under analysis. The result is judged according to the reference range of the reference group and the result beyond the normal range will be flagged.	NOTE If the Automatically match the customized reference group according to age and gender is set, gender and age of a patient will automatically match the reference group according to the corresponding relationship (No matter the reference group is selected or not). Refer to 5.4.3 Ref. Range for the setting of the reference group and range.
Department	Department receiving the patient.	Select from the dropdown list.
Area	Ward area of patient.	Select from the dropdown list or input directly.
Bed No.	Bed No. of inpatient.	Select from the dropdown list or input directly. NOTE The bed No. is required to be filled only for inpatients.

Parameter	Meaning	Operation
Sampling Time	Date and time when the sample is collected.	Click the date control for the settings. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd HH:mm, you should input the data in the sepuence of year, month, date, hour, and minute. Click or to select the date or
Time	sample is collected.	click the textbox to enter them directly. Click to delete the current data and re-enter information.
		NOTE The sampling time can be no later than the current
Submitter	Personnel submitting the sample.	system time. Select from the dropdown list or input directly.
Mode	Counting mode of the selected sample.	You do not need to enter it and it will be displayed automatically.
Delivery Time	Date and time when the sample is delivered.	 Click the date control for the settings. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd HH:mm, you should input the data in the sepuence of year, month, date, hour, and minute. Click or to select the date or click the textbox to enter them directly. Click to delete the current data and re-enter information. NOTE The delivery time can be no later than the current system time and cannot be earlier than the sampling time.
Checker	Personnel running the sample.	You do not need to enter it and it will be displayed automatically.
Run Time	Time when the sample is run.	You do not need to enter it and it will be displayed automatically.
Approver	Personnel validating the sample.	This parameter will be automatically displayed after the sample is validated.
Report Time	The date and time when the report is printed for the last time.	This parameter will be automatically displayed after the report is printed.

Parameter	Meaning	Operation
Diagnosis	Suspected diagnosis information.	Input in the textbox directly.
Remarks	Clarifications or notes.	Input in the textbox directly.

3. Click Apply or OK to save the changes.

8.4.8 **Query**

You can view the test results of a patient within a certain test date range by entering the query conditions. The procedures are shown as below:

1. Click the **Query** button to enter the multi-conditional query dialog box as shown below.

Query Sample ID Name Med Rec. No. 2015/07/23 2015/07/23 Run Date Sample status Not Validated Not Printed Not Transmitted Auto select All Samples OK Cancel

Figure 8-9 Query Conditions

2. Determine the query conditions as needed.

For the specific parameter description, see Table 8-2.

Table 8-2 Parameter Description of Query Conditions

Parameter	Meaning	Operation Description
Sample ID	Sample ID to be queried.	Input in the textbox.
Name	Name of patient.	Input in the textbox.
Med Rec. No.	Med Rec. No. of patient.	Input in the textbox.
Run Date	Test date range of sample.	Select the starting and ending dates of the sample test in the two data controls successively.

Parameter	Meaning	Operation Description
Sample status	Status of validation, printing or communication of the sample. Values: Not Validated Not Printed Not Transmitted	Please choose according to the actual situation. The default value is Not Validated .

NOTE

- Auto select checked by default indicates that the query result is being selected (with a blue background color). If it's unchecked, the query result will remain on a white background color.
- Click All Samples to close the current window, display all the samples again and restore all the filter conditions to the default values.
- 3. Click Query.

The system will display all the query results which meet the conditions.

8.4.9 Graph

In the **Review** interface, you can click **Graph** to browse the selected sample graph results, parameter results and flag messages. The procedures are shown as below:

- 1. Select a result to review in graph interface.
- 2. Click to unfold all function buttons.
- 3. Click **Graph** to enter the graph interface.

In the **Graph** interface, you can view sample information such as parameter results, graph results and flag messages. In addition, you can also print the analysis report as. See Figure 8-10.

Graph Next Print Previous Close Sample ID 1 Mode Venous Whole Blood Age Run Time 07-16-2015 11:27:32 Gender Name **WBC Message** Para. Result Unit Para. Result Unit **WBC** 6.57 10⁴3/uL RBC 4.62 10⁶/uL Lym% 37.4 % **HGB** 15.1 g/dL GRAN% 59.1 % **HCT** 46.7 % % MID% 3.5 MCV 101.1 fL **RBC Message** 10^3/uL 2.46 ↑ 32.6 Lym# MCH pg GRAN# 3.88 10⁴3/uL MCHC 32.3 g/dL 10^3/uL MID# 0.23 RDW-CV 13.3 % ↑ 57.5 RDW-SD fL **PLT Message** 10^3/uL PLT 284 MPV 9.0 fL PDW 15.5 PCT 0.256 % P-LCR 37.3 9/6 P-LCC ↑ 106 10⁹/L

Figure 8-10 Graphs Review

8.4.10 Microscopic Exam Para.

You can perform the microscopic examination parameter settings as per the following steps.

1. Click Microscopic Exam Para..

The microscopic examination parameters interface as shown in Figure 8-11 will pop up on the screen.

Microscopic Exam Para. Sample Type **Parameter Name** Value Capillary Neutrophilic segmented granulocyte 条 Exam. Time Neutrophilic band granulocyte 12-08-2015 Lymphocyte Microscopic Description Monocyte Eosinophil Basophil A Plasmacyte Atypical Lymph Blast Promyelocyte Neutrophilic myelocyte Eosinophilic myelocyte Basophilic myelocyte Apply Neutrophilic metamyelocyte Eosinophilic metamyelocyte $\overline{\mathbb{Y}}$ OK Basophilic metamyelocyte Cancel Prelymphocyte

Figure 8-11 Microscopic Exam Para.

2. Set the microscopic examination parameters by referring to Figure 8-3.

Table 8-3 Microscopic Exam. Parameters

Parameter	Meaning	Operation
Sample Type	Type of sample for microscopic examination. • Venous blood • Capillary	Click the Sample Type dropdown list box and select the type of sample for microscopic examination.
	 Cord blood 	
	Blood	

Parameter	Meaning	Operation	
Blood Type	Blood type of a patient.	Select the blood type of the patient in the Blood Type/ESR column. Click the first combo box next to the blood type, you can select from Blank , A , B , O and AB ; click the second combo box, you can select from Blank , RH+ and RH- .	
Exam. Time	Time of microscopic examination.	 Click the date control for the settings. The input sequence of the controls is the same with the date fomat on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd HH:mm, you should input the data in the sepuence of year, month, date, hour, and minute. Click or to select the date or click the textbox to enter them directly. Click to delete the current data and re-enter information. NOTE The Microscopic exam. time can be no later than the current system time. 	
Microscopic Description	Description of cells morphology.	Enter the morphology information for cells into the multi-line textbox.	

8.4.11 Customized Parameters

You can browse and edit the customized parameters results of the selected sample in the **Review** interface. The operation procedures are as shown below:

1. Click **Custom Para.** to enter the customized parameters setting interface as shown in Figure 8-12.




Figure 8-12 Customized Parameters

2. Click the cell corresponding to its Value column of the parameter, and enter the value.

If the unit and reference range of parameters have been set in the **Setup > Parameter > Custom Para.** interface, the corresponding unit and range (lower limit~upper limit) will be displayed in this tab. When both the value and range of parameters are numbers, and the number is out of the reference range, the relevant mark \uparrow or \downarrow will be displayed in the Flag column.

Please refer to **5.4.5 Customized Parameter** for customized parameters settings.

8.4.12 Run Chart

Operators can check and review run charts of sample parameter results in the database. There are three view modes: selected samples, samples on current page and samples on specified run dates.

- View the run chart of the selected sample (default)
 - a. Check no fewer than three sample records.
 - a. Click to unfold all function buttons.
 - b. Click Run Chart.

The system pops up a dialog box as shown below.

Run Chart

Selected Records

Samples on current page

Run Date

2015/08/14

OK

Cancel

Firgure 8-13 Viewing the Run Chart of the Selected Sample

c. Click OK.

The screen will show the parameter result run chart of the selected sample. See Figure 8-14.




Figure 8-14 Run Chart

- View the run chart of samples on current page
 - a. Click on the current page to unfold all function buttons.
 - b. Click the **Run Chart** button and select samples on current page in the pop-up dialog box. See Figure 8-15.

Run Chart

Selected Records

Samples on current page

Run Date

2015/07/23 - 2015/07/23

Figure 8-15 Viewing the Run Chart of Samples on the Current Page

c. Click OK.


The screen will show the parameter result run chart of the selected sample.

OK

Cancel

- View the run chart of samples on specified run dates
 - a. Click to unfold all function buttons.
 - b. Click the Run Chart button, and select Run Date in the pop-up dialog box.
 See Figure 8-16.

Figure 8-16 Viewing the Run Chart of Samples on Specified Run Dates



c. Click the date edit box, set a date range in the pop-up dialog box, then click **OK**.



- ➤ The input sequence of the controls is the same with the date format on the top right corner of the dialog box. For example, if the data format is yyyy/MM/dd, you should input the data in the sepuence of year, month, and date.
- Click or to select a date and time or enter the information in the textbox directly.
- > Click to clear the current data and re-enter the information.

d. Click OK.

The screen will show the parameter result run chart of the selected sample.

9 Quality Control

9.1 Introduction

Quality Control (QC) consists of strategies and procedures that measure the precision and stability of the analyzer. The results imply the reliability of the sample results. QC involves measuring materials with known, stable characteristics at frequent intervals.

Analysis of the results with statistical methods allows the inference that sample results are reliable. MR recommends running the QC program on a daily basis with low, normal and high level controls. A new lot of controls should be analyzed in parallel with the current lot prior to their Exp. dates. This may be accomplished by running the new lot of controls twice a day for five days using any empty QC file. The QC files calculate the mean, standard deviation and coefficient of variation for each selected parameter. The instrument-calculated results should be within the expected ranges published by the manufacturer.

NOTE

- You should only use the MR-specified controls and reagents. Store and use the controls and reagents by following the instructions for use of the controls and reagents. Controls beyond their Exp. date shall not be used. Controls (similar to standard blood samples) must be well mixed before use.
- General users only have the access for browsing and executing the QC analysis other than editing.

9.2 L-J Quality Control

9.2.1 QC Principle

In the L-J quality control, quality control can be applied to 21 parameters. You can set the QC information by setting the QC file before performing the QC analysis. Each QC file can be assigned 1 batch number for high, normal and low level controls. Each QC file can store up to 500 QC results. When there are more than 500 QC results, the new QC results will overwrite the oldest results in sequence.

9.2.2 QC Settings



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

NOTE

Only users with administrator-level access can edit the L-J settings.

Before running a new batch of controls, you need to assign a QC file to each batch of controls. You can complete the QC settings by setting QC information in the QC files.

9.2.2.1 Entering QC Information

The administrator can set the QC files by operations such as Copy, New, and Edit. Specific steps are shown below:

- 1. Click QC to access the QC interface.
- Click QC Settings to enter the QC Settings interface. See Figure 9-1.




Figure 9-1 L-J Quality Control

3. Click the New button or select a QC file (Existing/Total is 0/500) without QC counting results and

click the **Edit** button.

The interface as shown in Figure 9-2 will pop up on the screen.

File No. Limits (#) Limits (#) Para. Target Para. Target WBC PLT Lot No. MPV Lym% GRAN% PDW Level MID% PCT Normal v P-LCR Lym# Exp. Date GRAN# P-LCC 2015/07/23 MID# QC Mode RBC Whole Blood ₹ HGB QC Sample ID HCT MCV MCH Set Limits MCHC RDW-CV Save RDW-SD Close

Figure 9-2 Entering QC Information

You can also select the QC file of which data has been set and then click **Copy**, and edit the content based on the original data.

4. Set related information of the controls with reference to Table 9-1.

Table 9-1 QC File Information

Parameter	Parameter Description	Operation Description
File No.	QC file No Each QC file can store up to 500 QC results.	Read only.
Lot No.	Lot number of controls.	Enter into the textbox directly. NOTE The lot No. can not be empty and up to 16 digits can be entered. You can enter characters, numbers, letters and special characters, but no Chinese characters are allowed.
Level	Level of the controls, including 3 levels, i.e. High, Normal and Low	Select from the dropdown list.
Exp.Date	Exp. date of the controls.	The default Exp.Date is the current system date and needs to be changed to the actual Exp. date of the controls.

Parameter	Parameter Description	Operation Description
QC Mode	QC mode of the controls, including Whole Blood and Predilute.	Select from the dropdown list.
QC Sample ID	Users need to set the number of the controls here if he/she is used to performing the analysis with the controls placed among the daily samples. See section 9.2.3.2 Completing QC Analysis in the Sample Analysis Interface. If the user performs the analysis in the QC Analysis interface, the ID cannot be entered.	NOTE Letters, numbers and all characters that can be entered through the keyboard (including special characters) are allowed for the QC ID, but the number must end with a nonzero number. Chinese and other languages (such as Japanese, Korean, etc) are not supported. The length of the entries ranges from 1 to 25 and the entries shall not be empty. The last character of a sample ID must be numeric, but a string of "0" only is not an acceptable sample ID.
Target	Target of the QC parameter.	Enter the targets in the cell corresponding to the expected QC parameter according to the control target list with the corresponding lot No.
Limits (#)	Limits (#) of the QC parameter.	Enter the limits in the cell corresponding to the expected QC parameter according to the control target list with the corresponding lot No.

5. According to the target list of the corresponding lot No., enter the target and limits into the textboxes of the parameters to be included in the QC run.

You can set the display form of the limits or the calculation method of the limits among the preset values. Please refer to **9.2.2.2 Setting Limits** for the setting methods.


6. Click the Save button to save all the settings of the QC.

9.2.2.2 Setting Limits

You can take the following steps to adjust the display format of the limits and the calculation method of the preset limits.

- 1. Click **QC** to access the QC interface.
- 2. Click QC Settings to enter the QC Settings interface.
- 3. Click Set Limits.

A dialog box will pop up as shown below.




- 4. Select By SD or By CV according to the actual needs.
 - ➢ If By SD is selected, the limits will be displayed in form of absolute value.
 Click 2SD or 3SD to select either double or triple standard deviation to be the limits.
 - ➢ If By CV is selected, the limits will be displayed in form of percentage.
 Click the 2CV or 3CV to select either double or triple coefficient of variation to be the limits.
- 5. Click **OK** to save all the settings for the limits.

9.2.2.3 Deleting QC File

If you want to delete the QC files which will not be used any more, please take the following steps:

- 1. Click QC to access the QC interface.
- 2. Click QC Settings to enter the QC Settings interface.
- 3. Select the QC file to be deleted, and click **Delete**.

A dialog box will pop up as shown below.




4. Click Yes.

All selected QC files together with their QC results will be completely deleted.

9.2.2.4 Clearing QC results

If you want to delete QC results of a specified file, please take the following steps:

- 1. Click QC to access the QC interface.
- 2. Click QC Settings to enter the QC Settings interface.
- Select the QC file in which the QC results are expected to be deleted, and click Clear.A dialog box will pop up as shown below.



4. Click Yes.

QC results in the selected QC file will be deleted. As shown below. The value in the **Existing/Total** column will be restored to the initial value.

File No.	Lot No.	Level	Exp. Date	QC Mode	QC Sample ID	Editor	Existing / Total
1	L001	Normal	2015/12/30	Whole	QC001	admin	0/500
2	L002	Normal	2015/12/30	Whole	QC002	admin	0/500

9.2.3 Quality Control Analysis

After completing the QC settings, you can choose one of the following two modes according to the selected QC mode to run the quality control samples:

- Completing QC analysis in the QC Analysis interface
- Completing QC analysis in the Sample Analysis interface

9.2.3.1 Completing QC Analysis in the QC Analysis Interface



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.



WARNING

- The sample probe is sharp and potentially biohazardous. Exercise caution to avoid contact with the probe when working around it.
- The sample may spill from the unclosed collection tubes and cause biohazard. Exercise caution to the unclosed collection tubes.
- Collection tubes broken may cause personal injury and/or biohazard. Be sure to place the
 collection tubes in the right adapter before running, otherwise, the collection tubes may be
 broken and cause biohazard.
- Keep your clothes, hairs and hands away from the moving parts to avoid injury.
- The reagents are irritating to eyes, skin and mucosa. Wear proper personal protective
 equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when
 handling them in the laboratory.
- If the reagents accidentally spill on the skin, wash them off with plenty of water and if necessary, go see a doctor; if the reagents accidentally spill into the eyes, wash them off with plenty of water and immediately go see a doctor.



CAUTION

- Running quality controls in presence of errors may lead to incorrect analysis results. If you see
 the error alarms when running the quality controls, please stop and resume the analysis until the
 errors are removed.
- Do not re-use such disposable products as collection tubes, test tubes, capillary tubes, etc.
- Sample clump may lead to incorrect analysis results. Check if clump exists before running the controls; if it does, handle it as per the related laboratory procedures.

NOTE

- You should only use the MR-specified controls and reagents. Store and use the controls and reagents as instructed by instructions for use of the controls and reagents. Using other controls may lead to incorrect QC results.
- Before being used for analysis shake well the controls that have been settled for a while.
- Be sure to use the MR-specified disposable products including vacutainer blood collection tube, vacutainer blood collection tubes with anticoagulant and capillary tubes etc.

After completing the QC settings, users can perform the QC analysis in the **QC Analysis** interface Specific steps are shown below:

- 1. Click QC to access the QC interface.
- 2. Click **QC Analysis** and enter the QC analysis interface as shown in Figure 9-3.





Figure 9-3 QC Analysis

3. Select the QC file No. to be run.

The screen displays the corresponding file information, as shown in Figure 9-4.

Figure 9-4 QC File Information



- 4. Be sure that the level of the control to be run is the same with the current QC file, and the control to be run is not expired.
- Prepare the controls according to the set control mode and control instructions.
 Predilute the controls with reference to 6.5 Sample Collection and Handling and get diluted QC samples if the QC mode is Predilute.

NOTE

Be sure to evaluate predilute stability based on your laboratory's sample population and sample collection techniques or methods.

6. Shake the prepared control as shown below to mix it well.

Figure 9-5 Mixing the Controls



- 7. In the ready for counting state (namely, the indicator light of the main unit is green), place the controls under the sample probe where the probe can aspirate the well mixed sample.
- 8. Press the aspirate key and start running the controls.
- 9. Upon the completion of the aspiration, you'll hear a beep and you can remove the controls. When the running of QC analysis is complete, the QC results will be displayed in the current screen (as shown in Figure 9-6) and saved in the QC file automatically.

Sample Analysis Reagent Management Review QC 温息 Edit Restore QC Type L-J QC Analysis QC Graph QC Settings QC Table File No.: 1 QC Mode: Whole Blood Editor:develop Level:Normal Existing / Total:4/500 Exp. Date: 2015/12/31 Lot No.:L001 QC Sample ID:QC001 Run Time: 2015/08/14 18:22:43 Operator: admin Para. Result Unit Result Unit Para. WBC 1 6.57 10^3/uL PLT 284 10³/uL 37.4 9.0 0/0 fl 59.1 % 15.5 MID% 3.5 % 0.256 % 2.46 10⁴3/uL 37.3 % 3.88 10⁹/L 106 10⁹/L 0.23 10⁴9/L 100 200 RBC 4.62 10^6/uL PLT 15.1 g/dL 46.7 % fL 101.1 32.6 pq 32.3 g/dL RDW-CV RDW-SD 13.3 % fL 57.5

Figure 9-6 QC Analysis Results

10. Perform the above procedures to continue running the controls if necessary.

NOTE

- If the QC file is outdated, its valid period will be displayed in red.
- "↑" or "↓" alarm symbol will be displayed next to the results with deviations exceeding the set limits.

9.2.3.2 Completing QC Analysis in the Sample Analysis Interface



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.



WARNING

- The sample probe is sharp and potentially biohazardous. Exercise caution to avoid contact with the probe when working around it.
- The sample may spill from the unclosed collection tubes and cause biohazard. Exercise caution to the unclosed collection tubes.
- Collection tubes broken may cause personal injury and/or biohazard. Be sure to place the collection tubes in the right adapter before running, otherwise, the collection tubes may be broken and cause biohazard.
- Keep your clothes, hairs and hands away from the moving parts to avoid injury.
- The reagents are irritating to eyes, skin and mucosa. Wear proper personal protective
 equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when
 handling them in the laboratory.
- If the reagents accidentally spill on the skin, wash them off with plenty of water and if necessary, go see a doctor; if the reagents accidentally spill into the eyes, wash them off with plenty of water and immediately go see a doctor.



CAUTION

- Running quality controls in presence of errors may lead to incorrect analysis results. If you see
 the error alarms when running the quality controls, please stop and resume the analysis until the
 errors are removed.
- Do not re-use such disposable products as collection tubes, test tubes, capillary tubes, etc.
- Sample clumps may lead to incorrect analysis results. Check if clump exists before running the controls; if it does, handle it as per the related laboratory procedures.

NOTE

- You should only use the MR-specified controls and reagents. Store and use the controls and reagents as instructed by instructions for use of the controls and reagents. Using other controls may lead to incorrect QC results.
- Before being used for analysis shake well the controls that have been settled for a while.
- Be sure to use the MR-specified disposable products including vacutainer blood collection tube, vacutainer blood collection tubes with anticoagulant and capillary tubes etc.
- If the blood-sample mode is **Predilute**, then a reminder of predilute counting will pop up if the user presses the aspirate key to perform the counting. To close the prompt, please refer to **5.9 Auxiliary Settings.**

After completing the QC settings, you can place the controls among the daily samples and perform analysis together in the **Sample Analysis** interface. After the analysis is completed, the system will

store the results to the QC file with the corresponding ID.

Specific steps for performing QC analysis in the Sample Analysis interface are as follows:


Preparing the controls as per the set QC mode and instructions for the controls.
 Prediluting the controls with reference to 6.5 Sample Collection and Handling to get diluted QC samples if the QC mode is Predilute.

NOTE

Be sure to evaluate predilute stability based on your laboratory's sample population and sample collection techniques or methods.

2. Click Mode & ID in the Sample Analysis screen.

A dialog box will pop up as shown below.



- 3. Enter the set **QC Sample ID** in the **Sample ID** edit box (other options can be ignored). Refer to **9.2.2.1 Entering QC Information** for the setting of the **QC Sample ID**.
- 4. Well mix the prepared controls.
- 5. In the ready for counting state (namely, the indicator light of the main unit is green), place the controls under the sample probe where the probe can aspirate the well-mixed controls.
- 6. Press the aspirate key and start running the controls.
- Upon the completion of the aspiration, you'll hear a beep and you can remove the controls.
 When the running of the controls is complete, the QC results will be saved in the QC file automatically.
- 8. Perform the above procedures to continue running the controls if necessary.


NOTE

- If the QC file is outdated, its valid period will be displayed in red.
- "↑" or "↓" alarm symbol will be displayed next to the results with deviations exceeding the set limits.


9.2.3.3 Edit Result

Clicking **Edit** will allow you to edit the QC analysis result after the QC analysis is performed. See Figure 9-7.

Figure 9-7 Editing QC Results



The edited data will be marked with an **E**. As shown below.



9.2.3.4 Restore Result

Clicking **Restore** will allow the QC analysis results to be restored to the original results. After the data is restored, the **E** mark will disappear.

9.2.4 QC Result Review

After running controls, you can review the QC results in the following two forms:

- QC Graph
- QC Table

9.2.4.1 Graph



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

You can review the result of L-J QC graph as per the following steps.


- 1. Click **QC** to access the QC interface.
- 2. Click **Graph** and enter the interface as shown in Figure 9-8.

Reagent Management Sample Analysis Review QC Print Delete Outliers QC Type L-J QC Settings QC Analysis QC Graph QC Table File No.: 1 QC Mode: Whole Blood Level:Normal Existing / Total:3/500 Exp. Date: 2015/12/31 Editor:develop Lot No.:L001 QC Sample ID:QC001 Upper Limit Mean Para. 2015/08/14 18:22:39 Operator develop SD Target CV% Lower Limit 9.00 **WBC** 8.00 6.57 7.00 Lym% Lym# 111 MH

Figure 9-8 L-J QC Graph Interface

3. Select the QC file No. you want to review.

The screen will display the corresponding information and the graph. See Figure 9-9.



126

4. Click the buttons at the right side of the QC graph, then you can browse QC graphs of different parameters; click the buttons at the bottom of the QC graph, then you can browse all QC results.

Introduction to the Graph Interface

Upper Limit Operator admin 2015/08/17 11:30:44 Para. 3 SD Target Lower Limit 8.00 6.71 **WBC** 7.00 0.26 6.57 . 5 6.00 3.85 1 Lym% Lym# [44

Figure 9-10 L-J QC Graph Interface

Interface Description:

- 1 The Mean, SD and CV% of all the QC results of each parameter in the current graph.
- 2 The saving date and time of the QC points located on the gray line
- 3- The operator who run the QC analysis and obtained the QC points located on the gray line.
- 4 The QC results of the parameters that correspond to the QC points located on the gray line.
- 5 The QC points in each graph are displayed from left to right according to the sequence from the earliest to the latest. The QC points are connected by a line to illustrate the distribution trend.
- 6 The QC point corresponds to each QC result. Only the selected QC point displays its value under the parameter. The black QC point indicates the value is within the limit; the red QC point indicates the value is out of the limit.
- 7 When you clicking a QC point in the graph, the QC points of other parameters saved together with this one will be marked by a gray line.
- 8 The relative position of the QC point located on the gray line and the total QC points saved currently.

NOTE

The outliers are excluded from the calculation of Mean, SD and CV%.

Delete

The administrator can delete the QC results by the following steps:

- Delete a single QC result
 - a. Move the grayline to the desired QC result, and click **Delete**.
 - b. Select **Current Data** in the pop-up dialog box as shown in Figure 9-11.

Delete

Current Data

All Data

OK


Cancel

Figure 9-11 Deleting Current QC Data (QC Graph)

- c. Click OK.
- Deleting all the QC results in the current QC file

Click Delete, select All Data in the pop-up dialog box, then click OK. See Figure 9-12.

Figure 9-12 Deleting all QC Data (QC Graph)



Entering the Reasons for the Outliers

Do as follows to enter the reasons for the outliers:

1. Move the gray line to the desired QC point, and then click **Outliers**.

The pop-up window displays the QC results, reference values and deviation limits of all parameters corresponding to the gray line as shown in Figure 9-13.

The QC results exceeding the limit will be displayed in red.




Figure 9-13 Enter Cause of Outliers

- 2. You can select the reason from the given ones or manually enter the reasons (up to 200 characters) into the textbox after selecting Others.
- 3. Click **OK** to save the reasons for the outliers and exit.

NOTE

If you enter the reason for the group of QC points whose results are actually within the limits, then their corresponding QC data both in the QC Graph and QC Table will be displayed in red. And the data will return in black if you cancel the reason and then save the changes.

Print

You can have the QC data of the current page or all QC data in the QC file printed by clicking the **Print** button.

NOTE

The printed QC graph will not show any parameters which are not involved in the quality control.

9.2.4.2 Table



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

1. Click QC to access the QC interface.

2. Click **QC Table** to access the interface as shown in Figure 9-14.

Figure 9-14 L-J QC Graph Interface

	Export	Comm.	Print	ete	Dele			▼ L-J	L-J	QC Type
Ī			able	QC T	C Graph	sis Q	QC Analy	ettings	QC Se	
	al:4/500	xisting / Tota	E	lormal	Level:N	lood	de:Whole B	QC Mod		File No.: 3
	D:QC003	C Sample I	Q	:L003	Lot No.		admin	Editor: 8	2015/12/31	Exp. Date:
	MID%	MID#	GRAN%	GRAN#	Lym#	Lym%	WBC	Time	Date	
							7.00	I	Į.	Target
							1.00	1	1	Limits (#)
	↑ 3.5	↑ 0.23	↑ 59.1	↑ 3.88	† 2.46	↑ 37.4	6.57	11:30:44	2015/08	1
2	↑ 3.5	↑ 0.23	↑ 59.1	↑ 3.88	† 2.46	↑ 37.4	7.01	11:31:09	2015/08	2
	↑ 3.5	↑ 0.23	↑ 59.1	† 3.88	† 2.46	↑ 37.4	6.82	11:31:18	2015/08	3
	↑ 3.5	↑ 0.23	↑ 59.1	↑ 3.88	↑ 2.46	↑ 37.4	6.43	11:31:47	2015/08	4
Carre										


3. Select the QC file No. you want to review.

The screen will display the corresponding information and the graph.

4. Click the buttons at the bottom of the table to browse the QC data of desired parameters; Click the buttons on the right of the table to browse the QC results.

Editing

Choose a row in the QC table and click Edit, then you can edit the selected QC data.



The edited data will be marked with an **E**. See Figure 9-15.

Figure 9-15 Editing QC Results

	Date	Time	WBC
Target	1	1	7.00
Limits (#)	1	1	1.00
1	2015/08	11:30:44	6.57
2	2015/08	11:31:09	E 7.01

Restoring


Click **Restore** to cancel the editing of the QC results. After the data is restored, the **E** mark will disappear.

Delete

With the administrator-level access, users can delete the selected QC data, QC data on the current page and all QC data.


- Delete a selected QC result
 - a. Click the column containing the desired QC result, and then click **Delete**.
 - b. Select Current Data in the pop-up dialog box as shown in Figure 9-16.

Figure 9-16 Deleting Current QC Data (QC Graph)



- c. Click OK.
- Delete QC data on the current page
 - a. Click **Delete** on the page which contains the QC results expected to be deleted.
 - b. Select Current Data in the pop-up dialog box as shown in Figure 9-17.

Figure 9-17 Deleting all QC Data (QC Graph)




- c. Click OK.
- Delete all QC results

NOTE


Please note that this operation will delete all QC results of the selected QC file and cannot be reverted!

- a. Click Delete.
- b. Select All Data in the pop-up dialog box.



c. Click OK.

A dialog box will pop up as shown below.




d. Click Yes to delete all the QC results in the current QC file.

Print

All the QC data or the data within the specified date range can be transmitted to LIS/HIS. Specific steps are shown below:

- 1. Select a QC file No. to be printed.
- 2. Click Print.

A dialog box will pop up as shown below.




- 3. Select the QC data to be printed: all data or specified data.
 - Choose All Date to print all the QC date of the table.
 - Select Specified Data and set the date range in the date edit box, and QC data within the specified date range will be printed.
- 4. Click **OK** to print the data.

Export

If you wish to export the information and the result of the current QC file, do as follows:

- 1. Inset a USB flash disk in the USB interface on the machine.
- 2. Click Export.



A message box shown below will pop up.

3. Select the export path, and then enter the file name.

The file will be exported to the root directory of the USB flash disk (/udisk/sda1) and named in the format of SampleInfo_yyyyMMdd_hhmmss.csv. Among which, yyyyMMdd_hhmmss means data export year, month, date, hour, minute, and second.

4. Click Save.

When the export is finished, a message box as shown below will pop up.




Figure 9-18 Export successfully

5. Click **OK** to close the message box.

9.3 X-B Quality Control

9.3.1 QC Principle

The X-B analysis is a weighted moving average analysis that uses values obtained from patient samples. It uses the 3 red cell indices, MCV, MCH and MCHC to indicate the hematology instrument performance. This is QC with no controls, which is a method of performance control like QC with controls. Both methods reflect the analysis performance of the analyzer from different perspective. Thus, one method should not be replaced with the other.

It is recommended the X-B analysis be activated when the sample volume of your laboratory is greater then 100 samples per day. Effective use of X-B requires randomization of samples and a normal cross section of patients to prevent skewing of indices. A reference range is established by the given reference values as well as lower and upper limits for the purpose of observing the variation of QC results within the reference range.

The analyzer performs X-B QC for three parameters, MCV, MCH, and MCHC. Twenty to two hundred samples can be grouped together for X-B numerical analysis. The samples are derived from the results of normal analyzer counting, with no distinction of whole-blood or predilute mode. The analyzer can save maximum 1000 X-B QC results. When the saved QC results have reached the maximum number, the newest result will overwrite the oldest.

9.3.2 QC Settings



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

NOTE

Only users with administrator-level access can edit the L-J settings.

Perform the QC Settings before running the controls. You can complete the QC settings by entering the QC information.

9.3.2.1 Entering QC Information

You can complete the X-B QC settings as per the following steps:

- 1. Click **QC** to access the QC interface.
- 2. Select X-B from the dropdown list of the QC Type.
- 3. Click QC Settings.

You'll enter the QC Settings interface as shown in Figure 9-19.




Figure 9-19 X-B QC Setting

4. In the **Samples/Batch** edit box, enter the amount of samples to be included in calculating for an X-B QC point.

The range is between 20 and 200 and the recommended value is 20.

NOTE

Once the **Samples/Batch** is changed, the number of valid sample results will be recalculated. For example, if 20 valid samples are needed for the X-B QC calculation, when you change the value of Samples/Batch after 10 group of valid sample results have been acquired, these 10 group of results will be discarded, and only valid sample results generated afterwards will be used in the QC calculation.

- 5. Click the **Open** button of **X-B** to open the X-B quality control. The samples results will be included to calculate the X-B.
- 6. Enter the targets and limits for the QC parameters.

NOTE

- All the targets and limits for the QC parameters must be entered.
- When first use, the default setting will provide the Initial values for the targets and limits of the three QC parameters.
- If the QC data have existed in the QC file, you are not allowed to edit the target and limits.

You can set the display form of the limits or the calculation method of the limits among the preset values. See section **9.3.2.2 Setting Limits.**

7. Set the valid upper and lower limits for the QC parameter in **Sample Validity Setting** field. Setting sample validity is to set the valid range of four QC parameters, RBC, MCV, MCH and MCHC. To be incorporated into X-B QC calculation, the sample results should satisfy the validity ranges of all these four parameters.

NOTE

Once the **Samples/Batch** is changed, the number of valid sample results will be recalculated. For example, if 20 valid samples are needed for the X-B QC calculation, when you change the value of Samples/Batch after 10 group of valid sample results have been acquired, these 10 group of results will be discarded, and only valid sample results generated afterwards will be used in the QC calculation.

8. Click the **Save** button to save all the settings of the QC.


If the entered value exceeds the acceptable range or the upper limit is lower than the lower limit, a reminder message will pop up and you will be prompted to re-enter the correct data and save the entry again.

9.3.2.2 Setting Limits

You can take the following steps to adjust the display format of the limits and the calculation method of the preset limits.

1. Click Set Limits.

A dialog box will pop up as shown below.



- 2. Select By SD or By CV according to the actual needs.
 - ➢ If By SD is selected, the limits will be displayed in form of absolute value.
 Click 2SD or 3SD to select either double or triple standard deviation to be the limits.
 - ➢ If By CV is selected, the limits will be displayed in form of percentage.
 Click the 2CV or 3CV to select either double or triple coefficient of variation to be the limits.

3. Click **OK** to save all the settings for the limits.

9.3.2.3 Restoring Defaults

In QC setting, click **Restores Defaults** button to restore the parameter reference values, limits and sample validity to the default settings.

NOTE

- If the QC data have existed in the QC file, you are not allowed to restore the parameters.
- Clicking Restores Defaults can only store the default settings of Target, Limits and Sample Validity Setting, while Samples/Group, X-B QC switch and limit settings cannot be restored.

9.3.3 Quality Control Analysis



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

After the QC settings, the analyzer will automatically start the X-B QC analysis.

After every 20~200 results (determined by the setting) are obtained, the system will perform the X-B calculation once automatically. You can review the result in X-B graph or X-B table.

In X-B QC, sample results conforming to any of the following conditions will be considered as invalid and can not be used in the QC calculation.

- Sample results exceeding the linearity range
- Background results
- Sample results not conforming to the Sample Validity Setting
- QC data for other QC programs (such as L-J QC)
- Calibration data
- Results generated while there are errors which could affect the accuracy of the results (insufficient aspiration volume or clogging for example).

9.3.4 QC Result Review

After running controls, you can review the QC results in the following two forms:

- QC Graph
- QC Table

9.3.4.1 Graph



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

Access the X-B QC Graph interface by taking the following steps:

- 1. Click **QC** to access the QC interface.
- 2. Select X-B from the dropdown list of the QC Type.
- 3. Click Graph.

The X-B QC Graph interface will be displayed. See Figure 9-20.





Figure 9-20 QC Graph

4. You can also drag the scroll bar down to the graph horizontally to browse all the QC results.

Introduction to the Graph Interface

Figure 9-21 X-B QC Graph




- 1 The amount of samples included in calculating for each QC point.
- 2 The saving date and time of the QC points located on the green line
- 3 The Mean, SD and CV% of all the QC results of each parameter in the current graph.
- 4 The QC results of the parameters that correspond to the QC points located on the green line.
- 5 The QC points in each graph are displayed from left to right according to the sequence from the earliest to the latest. The QC points are connected by a line to illustrate the distribution trend.
- 6 The QC point corresponds to each QC result. Only the selected QC point displays its value under the parameter. The black QC point indicates the value is within the limit; the red QC point indicates the value is out of the limit.
- 7 When you clicking a QC point in the graph, the QC points of other parameters saved together with this one will be marked by a green line.
- 8 The relative position of the QC point located on the green line and the total QC points saved currently.

Delete

The administrator can delete the QC results by the following steps:

- Delete a single QC result
 - a. Move the green line to the desired QC result, and click Delete.
 - b. Select Current Data in the pop-up dialog box as shown in Figure 9-22.


Figure 9-22 Deleting Current QC Data (QC Graph)



- c. Click OK.
- Deleting all the QC results in the current QC file

Click **Delete**, select **All Data** in the pop-up dialog box, and then click **OK**. See Figure 9-23.

Figure 9-23 Deleting all QC Data (QC Graph)



Print

Click the **Print** button to print the QC table.

9.3.4.2 Table



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

Access the X-B QC Graph interface by taking the following steps:


- 1. Click QC to access the QC interface.
- 2. Select X-B from the dropdown list of the **QC Type**.
- 3. Click QC Table.

The X-B QC table interface will be displayed. See Figure 9-24.

Sample Analysis Reagent Management Review QC Print Comm. Export Delete QC Type X-B QC Table QC Settings QC Graph Editor Samples/Group 1/Group admin Date Time MCV MCH MCHC Target 1 1 89.5 30.5 34.0 1 1 Limits (#) 2.7 0.9 1.0 2015/08/14 18:06:04 1 91.1 30.6 34.5 30.9 2 2015/08/17 15:50:55 101.1 33.3

Figure 9-24 QC Table

Introduction to the QC Table Interface




- 1 The amount of samples included in calculating for each QC point.
- 2 QC parameters (displayed in the same order as the **Graph** screen)
- 3 The No. of the QC result saved in the QC file (arranged from left to right in the order that from the earliest to the latest)
- 4 QC Result. The value of the QC result is the X-B result of each batch of samples.
- 5 QC flag: The flag ↑ or ↓ will be used to prompt the result that out of the limits

Delete

With the administrator-level access, users can delete the selected QC data, QC data on the current page and all QC data.


- Delete the selected QC data
 - a. Choose the desired QC result, and then click **Delete**.
 - b. Select Current Data in the pop-up dialog box as shown in Figure 9-25.

Figure 9-25 Deleting Current QC Data (QC Graph)



- c. Click OK.
- Delete current page QC data
 - a. Click **Delete** on the page which contains the QC results expected to be deleted.
 - b. Select **Current Page Data** in the pop-up dialog box as shown in Figure 9-26.

Figure 9-26 Deleting all QC Data (QC Graph)




c. Click OK.

Delete all QC data

NOTE


Please be careful to perform this operation as it will delete all QC data of the selected QC file.

- a. Click Delete.
- b. Select All Data in the pop-up dialog box.



c. Click OK.

A dialog box will pop up as shown below.




d. Click Yes to delete all data of the current QC file.

Print

You can print all the QC data or the data within the specified date range of the selected QC file. Specific steps are shown below:

- 1. Select a QC File No. to be printed.
- 2. Click Print.

A dialog box will pop up as shown below.




- 3. Select the QC data to be printed: all data or specified data.
 - > Select All Data to print the data in the whole QC list.
 - > Select **Specified Data**, set the date range in the date edit box, and QC data within the specified date range will be printed.
- 4. Click **OK** to print.

Export

If you wish to export the information and the result of the current QC file, do as follows:

- 1. Insert a USB flash disk in the USB interface on the analyzer.
- 2. Click Export.

A message box shown below will pop up.




3. Select an export path for the data and enter the file name.

The file will be exported to the root directory of the USB flash disk (/udisk/sda1) and named in the format of SampleInfo_yyyyMMdd_hhmmss.csv. Among which, yyyyMMdd_hhmmss means data export year, month, date, hour, minute, and second.

4. Click Save.

When the export is finished, a message box as shown below will pop up.

Figure 9-27 Export successfully



5. Click **OK** to close the message box.

10 Calibration

10.1 Introduction

Calibration is a procedure to standardize the analyzer by determining its deviation, if any, from calibration references and to apply any necessary correction factors. To get accurate blood analysis results, perform calibration of the analyzer following the procedures given in this chapter when it's needed.

NOTE

- Calibration procedures can only be performed by users with the administrator-level access. The login users with the access level of general users can not perform the calibration procedures but only browse the calibration coefficients.
- You should only use the MR- specified calibrators and reagents. Store and use the calibrator and reagents following the instructions for use of the calibrations and reagents.
- The analyzer identifies a sample as a calibration sample only if the analysis is started from the Cal interface.
- The calculation of repeatability is included in the calibration procedure.

10.2 When to Calibrate

This analyzer is calibrated at the factory just before shipment. It is electronically stable and does not require frequent recalibration if you operate and maintain it as instructed by this manual. You need to recalibrate this analyzer if:

- it is the first time this analyzer has been used (usually done by a MR-authorized representative when installing the analyzer).
- an analytical component has been changed.
- the quality control results indicate that there may be a problem.
- the operating environment (such as the temperature) has changed significantly.

NOTE

- All of the measured parameters must be calibrated before readings of this analyzer can be used as valid analysis results.
- For laboratories conducting routine tests, the calibration should be applied at least once every six months.

10.3 How to Calibrate

There are three calibration programs available on this analyzer: manual calibration, auto calibration

using calibrators and auto calibration using fresh blood samples.

All or part of the parameters of WBC, RBC, HGB, MCV and PLT can be calibrated by the calibration procedure.

10.3.1 Preparation



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.



WARNING

- The sample probe tip is sharp and may contain biohazardous materials. Exercise caution to avoid contact with the probe when working around it.
- The reagents are irritating to eyes, skin and mucosa. Wear proper personal protective
 equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when
 handling them in the laboratory.
- If the reagents accidentally spill on the skin, wash them off with plenty of water and if necessary, go see a doctor; if the reagents accidentally spill into the eyes, wash them off with plenty of water and immediately go see a doctor.
- Keep your clothes, hairs and hands away from the moving parts to avoid injury.
- Be sure to dispose of reagents, waste, samples, consumables, etc. according to local legislations and regulations.



CAUTION

Do not re-use such disposable products as collection tubes, test tubes, capillary tubes, etc.

NOTE

- You should only use the MR-specified controls and reagents. Store and use the controls and reagents as instructed by the instructions for use of the controls and reagents.
- Be sure to use the MR-specified disposable products including vacutainer blood collection tube, vacutainer blood collection tubes with anticoagulant and capillary tubes etc.

Carry out the calibration only when the background range, repeatability and carryover are within the specified limits given in the manual, otherwise, the problems must be identified and solved before you determine if calibration is needed. If you cannot solve the problems, please contact MR Service Department.

- 1. Check and make sure enough reagents have been prepared for the calibration. You need to start over the calibration if the reagents run out during the process.
- 2. Do the background check.
 - If the analyzer alarms are activated for **abnormal background** results, see **13 Troubleshooting** for solutions. (Refer **A.4.2 Normal Background** to for background range.)
- 3. Run the median controls in Whole Blood mode consecutively for 11 times, take and view repeatability of the counting results from the 2nd run through the 11th run in the **Review** interface

and make sure they are within the range specified in A.4.4 Repeatability.

4. Run the corresponding diluent for 3 times immediately after running the high-level controls for 3 times and calculate the carryover by the following formulae:

The calculated carryovers shall meet the requirements in A.4.5 Carryover.

5. It is recommended that you create a log table for your analyzer. This log table should contain all the necessary information pertinent to your analyzer. The suggested items that you may want to include in the log table are: calibration date, supplier of calibrator, lot number, expected results and limits, and result of background check.

10.3.2 Manual Calibration

Complete the manual calibration as per the following procedure:

- 1. Click **Cal** in the menu page to access the calibration interface.
- 2. Click **Manual** to access the manual calibration interface. See Figure 10-1.




Figure 10-1 Manual Calibration

The calibration coefficients of whole blood mode and predilute mode are displayed on the **Manual** interface.

NOTE

The login users with the access level of general users can not perform the calibration procedures but only browse the calibration coefficients on the current screen. To perform the calibration, please log out and then log in as users with administrator-level access.

3. Check the calibration coefficient and calculate the new coefficient using the following equation.

$$\begin{tabular}{ll} New calibration factor & \underline{ Current calibration factor \times \\ Mean & \underline{ Ref erence value } \\ \hline \end{tabular}$$

For example, the WBC reference value of a calibrator is 8.3, and the current calibration coefficient of the whole blood mode is 99.00%.

Run the calibrator in whole blood mode for 11 consecutive times and calculate the WBC results of the 2nd to 11th runs (n=10): 8.4, 8.2, 8.2, 8.3, 8.3, 8.1, 8.2, 8.1, 8.2, 8.2. The obtained CV is 1.1% and the Mean is 8.22, which meet the requirements.

The new calibration coefficient is obtained:

New calibration factor=
$$\frac{99.00\% \times 8.3}{8.22}$$
=99.96%

The calculated calibration coefficients shall be between 75%~125%. In case of an invalid calibration coefficient, try to find out the reason (e.g. calibration material not thoroughly mixed, incorrect operation, etc.). Then recalibrate the analyzer and recalculate the calibration coefficients.

4. Enter the new calibration coefficients into the factor cell of the parameter that requires calibration.


NOTE

The entered calibration coefficients shall be between 75.0%~125.0% (calculation results rounded to two decimal places).

5. Click Save.

➤ If the new calibration coefficient is valid and different from the original value, the following dialog box will pop up.3


Figure 10-2 Calibration set successfully



On the screen, the calibration coefficient is refreshed to be the new one and the calibration date is refreshed to be the current system date.

If the new calibration coefficients are invalid, the message box will pop up. Click **OK** to close the message box and enter a valid factor.

Figure 10-3 Invalid Coefficients



- 6. (Optional) Click Print to print the current calibration coefficient.
- 7. Click **Exit** to close the **Manual** interface.

10.3.3 Auto Calibration Using Calibrators



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

NOTE

- Only MR-specified calibrators shall be used. MR will not be responsible for any erroneous result caused by using other calibrators.
- See the instructions for use of the calibrators for the lot No., Exp. Date and the target.

Complete the calibration with calibrators as per the following procedure:

- 1. Click **Cal** in the menu page to access the **calibration** interface.
- 2. Click Calibrator.

The Calibrator interface pops up as shown in Figure 10-4.




Figure 10-4 Auto Calibration Using Calibrators

- 3. Enter the lot No. of the calibrator into the Lot No. box.
- 4. Click the **Exp. Date** box, and then edit the Exp. Date.

NOTE

- The Exp. Date can be no earlier than the current system date.
- The entered Exp. Date should be either the Exp. Date printed on the labeling or the open-container Exp. Date, whichever is earlier. The open-container Exp. Date is calculated as follows: the date on which the container is opened + the open-container stability days.
- 5. Input the target values of the parameters in the corresponding cell.
- 6. Prepare the calibrators following their instructions for use and place the calibrators under the sampling probe.
- 7. Press the aspirate key to start the calibration counting.
 - > The valid results within the linearity range will be displayed directly.
 - If the calibration counting data of any parameter in the current counting are out of the display range or linearity range of the parameter, a message box will pop up on the screen prompting that the calibration data is invalid.
 - Click **OK** to close the message box and delete the data from the table without saving.
 - ➤ If any of the parameter's value in the calibration counting differs from the Target value by more than 50%, the system will prompt you with a message box asking if the calibration counting results should be kept.
 - To keep the results, click Yes; to remove the results, click No.

NOTE

- After the valid calibration result is obtained, the parameters with corresponding checkboxes ticked off will be involved in the calculation of the calibration coefficients by default.
- If you switch to other interfaces before the new calibration coefficients are obtained, the system will discard the current calibration data and keep the original calibration coefficients.
- 8. To get 10 valid counting results, repeat steps 6~7 ten times.
 - The analyzer will, by default, calculate the Mean, CV% and the new calibration coefficients based on all the ticked-off calibration data according to the formulae.
- 9. You can select a few groups of data for the calculation of the calibration coefficients which can be obtained unless at least 5 groups of ticked-off data are included. Each time when you tick off or uncheck the checkboxes, the calibration coefficients will be refreshed and displayed in time.
 - When the amount of the valid calibration data in the list reaches 10, a message box of **Calibrator calibration done!** will pop up. Click **OK** to close the message box.

NOTE

The out-of-range CV% does not influence the display of the calibration coefficients.

10. Click Save.

- ➢ If the calculated calibration coefficient is within the range of 75%~125% (i.e.,>=75% and <=125%) and the CV% values of all the calibration parameters are not beyond the repeatability index, a dialog box prompting the successful calibration setting will pop up. Click OK to close the message box.</p>
- ➢ If the obtained calibration coefficient of any parameter is not within the range of 75%~125% or the CV% of any calibrated parameter does not meet the repeatability, the calibration coefficient will not be saved and a dialog box will pop up. Click Yes to close the dialog box and repeat the calibration operations.
- 11. (Optional) Click **Print** to print the calibration results.

10.3.4 Auto Calibration Using Fresh Blood Samples



All the samples, controls, calibrators, reagents, wastes and areas in contact with them are potentially biohazardous. Wear proper personal protective equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when handling relevant items and areas in the laboratory.

Complete the calibration using fresh blood samples as per the following procedure:

- 1. Click **Cal** in the menu page to access the calibration interface.
- 2. Click Fresh Blood.

The fresh blood sample calibration interface pops up, as shown in Figure 10-5.




Figure 10-5 Auto Calibration Using Fresh Blood Samples

- 3. Prepare 3 to 5 normal fresh blood samples as instructed by **6.5 Sample Collection and Handling**.
- 4. Run each of the prepared samples on the reference instrument three times at least. Average the results for your reference values.

NOTE

The reference instrument must be a properly running standard analyzer so as to ensure the accuracy of the reference values.

- 5. Enter the reference values for the parameters to be calibrated in the corresponding **Target** textbox.
- 6. Place the blood sample under the sampling probe, press the **aspirate key** on the analyzer to run the samples.

The system will calculate the values for WBC, RBC, HGB, MCV and PLT of the sample.

7. Repeat step 6 for 10 times and calculate the counting results for sample No. 1 in the 10 runs. The system will calculate the Mean, CV and Calibration coefficient for each parameter of the sample.

If the obtained calibration coefficient for any sample is not within the valid range or CV% or any calibrated parameters does not meet the repeatability, a dialog box indicating invalid new calibration coefficient will pop up when you are selecting other blood samples.

- Click **Yes** to clear the calibration data of the sample. Redo the calibration or redo after running another sample meeting all criteria.
- 8. Refer to steps 6~7 and perform the counting operations for the remaining four blood samples. The system will calculate the Mean, CV and Calibration Coefficient for each parameter of the remaining 4 blood samples.
- 9. Click Calculate.

The system will calculate the average of the calibration coefficients, namely, the mean calibration coefficient (%), as the new calibration coefficient based on the five blood samples.

You can also check at least three accurate calibration coefficients and the system will re-calculate the mean calibration coefficient (%).

NOTE

The mean calibration coefficient is invalid if its absolute value of deviation from the original calibration coefficient is greater than or equal to 5%.

10. Click Save.

- ➢ If the mean calibration coefficient is within the valid range (the absolute value of deviation from the original calibration coefficient is less than 5%), you'll be prompted that the mean calibration coefficient is saved successfully.
- If the mean calibration coefficient is not within the valid range (the absolute value of deviation from the original calibration coefficient is greater than or equal to 5%), you'll be prompted that the mean calibration coefficient is invalid.

NOTE

CV% out of standard will not affect the display of calibration coefficient.

- 11. Click **OK** to close the message box.
- 12. (Optional) Click **Print** to print the calibration results.

10.4 Verifying Calibration Coefficients

It is recommended that you take the following steps to verify the calibration coefficients:

- 1. Run the calibrator at least three times and check whether the means of the obtained results are within the expected ranges.
- 2. Run the low-, normal- and high-level controls each for three times at least, and check whether the means of the obtained results are within the expected ranges.
- 3. Run at least three fresh blood samples with known reference values, each for six times at least, and check whether the means of the obtained results are within the expected ranges.

11 Reagent Management

Once the new reagent is connected to the analyzer, you should set the reagent configurations, including validity period, residue volume and reagent barcode on the **Reagent Management** interface. Upon the completion of reagent configuration, you can perform the procedures for reagent replacement.



WARNING

- The reagents are irritating to eyes, skin and mucosa. Wear proper personal protective
 equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when
 handling them in the laboratory.
- If the reagents accidentally spill on the skin, wash them off with plenty of water and if necessary, go see a doctor; if the reagents accidentally spill into the eyes, wash them off with plenty of water and immediately go see a doctor.


NOTE

- After long-distance transportation, the reagent must be allowed to settle for more than one day before use.
- When you have replaced the diluent, cleansers or lyses, run a background check to see if the results meet the requirement.

11.1 Accessing the Interface

Click **Reagent Management** in the menu navigation area or menu page, to access the reagent management setting interface. See Figure 11-1.

Figure 11-1 Reagent Management



Refer to Table 11-1 for related parameter descriptions.

Table 11-1 Parameter Description for Reagent Management

Parameter	NOTE	
	Current model of the analyzer.	
	Open system	
Current Model	Closed system	
	Reagent setting procedures for different analyzer models vary, please refer to 11.2 Reagent Information Settings .	
Reagent Name	Name of the reagent.	
Exp.Date	Exp. Date of the unopened reagent will be shown upon the completion of the reagent settings. Any reagent, regardless of its container being opened or not,	
	should not be used beyond this date.	
Open-container Date	The date on which the reagent container is opened. The default open-container date is the date on which the reagent settings are completed.	
Period after opening (PAO)	The validity period (days) after the reagent container is opened. It will be shown upon the completion of the reagent settings.	
Open-container Exp. Date	Exp. Date of the opened reagent, and it will be shown upon the completion of the reagent settings.	
Residue Volume	The current residue volume of the reagent, and it will be shown in ml upon the completion of the reagent settings. The unit is ml.	

11.2 Reagent Information Settings

Once the new reagent is connected to the analyzer, you need to set the reagent configurations, including validity period, residue volume and reagent barcode on the **Reagent Management** interface. Upon the completion of reagent configuration, you can perform the procedures for reagent replacement.

Reagent setting procedures for different analyzer models vary. The reagent setting procedures for both open and closed models will be presented on the following pages.

11.2.1 Open system

For open systems, reagent setting procedures are as follows:

1. Select the reagent to be set, and then click **Setup**.

This launches the **Reagent Information** page as shown in Figure 11-2.

Reagent Information
Reagent Name
DIL-E Diluent
Exp. Date
//
Period After Opening

Residue Volume

L

DIL-E Diluent
Exp. Date
//
Period After Opening

Barcode 1
Barcode 2

Apply

Figure 11-2 Reagent Information

- 2. To enter the reagent information, use any of the following methods.
 - Manual Entry

Detailed parameter description is shown in Table 11-2.

Table 11-2 Parameter Description of Reagent Information

Table 11-2 Parameter Description of Reagent Information				
Parameter	Meaning	Operation		
Reagent Name	Name of the reagent to be set.	Input in the textbox directly.		
Exp.Date	The expiration date of the unopened reagent (see the outer packaging of the reagent). Any reagent, regardless of its container being opened or not, should not be used beyond this date.	 Click the date control for the settings. The input sequence of the controls is: year, month, and date. Click or to select the date or click the textbox to enter them directly. Click to delete the current data and re-enter information. NOTE The validity date of the reagent can be no later than the validity date indicated on the packaging and cannot be earlier than the current system date. 		
Period after opening (PAO)	The validity period (days) of the open-container reagent (see the product packaging).	Input in the textbox directly.		

Parameter	Meaning	Operation
Residue Volume	The current residue volume of the reagent (ml).	Input in the textbox directly.

Manually input the reagent barcode, and click **Load**; or input the barcode via a peripheral barcode scanner.

A correctly entered barcode will prompt a message shown below the barcode box, indicating a successful loading, and the validity date and residue volume will be shown in the corresponding textboxes.

If the bar code fails to be loaded, check if the reagent has been used or expired and the reagent name is correct. If all the information is correct, but the failure persists, please contact MR After-sales Service Department.


3. Click Apply.

The system message will pop up, indicating the successful reagent settings.

DIL-E Diluent set successfully!

Figure 11-3 Successful Reagent Settings

- 4. Click OK.
- 5. Continue to perform 1~4 and set the other reagent information; or click interface.



NOTE

- Once the reagent settings are successfully completed, the system prompt at the top right corner
 of the screen will show that the reagent has not been replaced. To remove this error, click the
 error message and then click **Remove Error** in the pop-up dialog box. The analyzer will
 complete the reagent replacement and remove the error automatically.
- When you have replaced the diluent, cleansers or lyses, run a background check to see if the results meet the requirement.

11.2.2 Closed system

For closed systems, reagent setting procedures are as follows:

1. Select the reagent to be set, and then click **Setup**.

This launches the Reagent Information page as shown in Figure 11-4.




Figure 11-4 Reagent Information

2. Input the barcode via a peripheral barcode scanner or manual input, and then click **Load**.

A correctly entered barcode will prompt a message shown below the barcode box, indicating a successful load, and the validity date and residue volume will be shown in the corresponding textboxes.

If the barcode fails to be loaded, check if the reagent has been used or expired and the reagent name is correct. If all the information is correct, but the failure persists, please contact MR After-sales Service Department.

3. Click Apply.

➤ For the settings of diluents, a pop-up dialog box as shown in Figure 11-5 indicates the completion of the settings. Please perform steps 5~6.





Figure 11-5 Successful Diluent Settings

For the settings of lyses, a dialog box as shown in Figure 11-6 will pop up. Please perform the next step.


Figure 11-6 IC Card Verification



4. Put the RF card attached to reagent packing on the RF card reader in front of the analyzer.

The beeping of the card reader and a pop-up dialog box as shown in Figure 11-7 indicate the successful reagent settings.

Figure 11-7 Successful Reagent Settings



NOTE

- The IC card is intended for single use only.
- If IC card verification fails, please follow the system prompts and use a valid IC card for re-reading.
- 5. Click OK.
- 6. Click Close to exit.

NOTE


- After the reagent settings are successfully completed, the error message area at the top right of
 the software interface will prompt that the reagent is not replaced. Please click the error
 message and click Remove Error in the pop-up dialog box. The analyzer will complete the
 replacement of the reagent and remove the error.
- When you have changed the reagents, run a background check to see if the results meet the requirement.

11.3 Reagent Replacement

After completing the reagent settings, you should perform the reagent replacement operations. You can select to replace one type of reagent at a time or all reagents. The statistical method is applied as follows:

1. Select a type of reagent to be replaced, and click **Replace**; or click **Replace All** to replace all the reagents.

After the replacement is completed, a message box as shown below will pop up on the screen.



2. Click **OK** to close the message box.

12 Service

12.1 Introduction

This analyzer provides multiple maintenance functions for this purpose. This chapter introduces how to use the provided functions to maintain and troubleshoot your analyzer. Preventive and corrective maintenance procedures are required to keep the analyzer in a good operating condition.



All the analyzer components and surfaces are potentially infectious, take proper protective measures for operation or maintenance.



CAUTION

- Performing unauthorized maintenance procedures can damage your analyzer. Do not perform any maintenance procedures that are not described in this chapter.
- In case of problems not specified in this manual, contact MR customer service department or your local agent for assistance.
- Only MR-supplied parts can be used for maintenance. For any question, contact MR customer service department or your local agent.
- Exercise caution to avoid contact with the sharp sample probe when performing maintenance.

12.2 Maintenance

The analyzer provides multiple service functions helping users to perform daily maintenance.

12.2.1 Reagent Replacement



WARNING

- The reagents are irritating to eyes, skin and mucosa. Wear proper personal protective
 equipment (e.g. gloves, lab uniforms, etc.) and follow laboratory safety procedures when
 handling them in the laboratory.
- If the reagents accidentally spill on the skin, wash them off with plenty of water and if necessary, go see a doctor; if the reagents accidentally spill into the eyes, wash them off with plenty of water and immediately go see a doctor.

NOTE

- After long-distance transportation, the reagent must be allowed to settle for more than one day before use.
- When you have replaced the diluent, cleansers or lyses, run a background check to see if the results meet the requirement.

You should replace the reagents when:

- The system indicates that the reagent is used up
- The suspicious flag indicates that the reagent in the pipeline is contaminated
- The reagent is contaminated or expired
- WBC or RBC bubbles are identified.


You can replace any of the following reagents:

- DIL-E Diluent
- LYE-1 LYSE

Do as follows to replace the reagents:

- 1. Refer to Figure 2-2 in **2.6.1 Electrical Connections** for reagent connections.
- 2. Click the **Service** icon in the menu page to access the **Service** interface as shown in Figure 12-1.

Figure 12-1 Service



3. Click Replace Reagent in the Maintenance selection.

The interface as shown in Figure 12-2 will pop up on the screen.

Figure 12-2 Reagent Replacement



Double click the name of the reagent that needs to be replaced, such as Replace DIL-E Diluent.
 After the replacement is completed, the following message box will pop up.

Figure 12-3 Reagent Replaced



- 5. Click **OK** to close the message box.
- 6. Perform the above procedures to replace other reagents if necessary.

12.2.2 Cleaning

Clean corresponding parts according to the actual situation:

WBC bath

When the background of WBC- and/or HGB-specific parameters exceeds the Ref. Range, you should clean the WBC bath.

RBC bath

When the background of RBC- and (or) PLT-specific parameters exceeds the Ref. Range, you should clean the RBC bath.

Sample probe

When the sample probe is dirty, you should clean the sample probe.

The cleaning procedures are as follows:

1. Click the **Service** icon in the menu page to access the **Service** interface.




Figure 12-4 Service

2. Click **Clean** in the **Maintenance** selection, an interface as shown in Figure 12-5 will pop up on the screen




Figure 12-5 Cleaning

3. Double click the icon of the part that needs to be cleaned, such as Clean Sample Probe.

When the system cleaning is complete, the message box will pop up to show that the cleaning is done.

Figure 12-6 Cleaning Done



- 4. Click **OK** to close the message box.
- 5. Perform the above procedures to clean other components if necessary.

12.2.3 Maintenance

Maintenance of the analyzer includes unclogg and cleanser soak.

12.2.3.1 Unclogging

If clogging is found, or it is suspected that the counting results are not accurate due to aperture clogging, you can perform the unclogging operations.

The unclogging procedures are shown as follows:

1. Click the **Service** icon in the menu page to access the **Service** interface.

Service Self-test Maintenance Status Log Replace Reagent Syringe Self-test All Logs Temperature Clean Pressure Self-test Voltage/Current Set Paras Other Self-test Maintain Version Info. Fault Logs Comprehensive Valve/Pump Self-test Other Logs Device Cal Touch Screen Cal.

Figure 12-7 Service

2. Click Maintain in the Maintenance selection.


The interface as shown in Figure 12-8 will pop up on the screen.

Figure 12-8 Maintenance



3. Click the Unclog icon.

The system will start clogging, and a message box will pop up. After the unclogging is completed, a message box will pop up to show that the clogging is done.



- 4. Click **OK** to close the message box.
- 5. Perform the above procedures to continue unclogging if necessary.

12.2.3.2 Cleanser Soak


The cleanser soak should be performed under the following circumstances:

- When the problems including the background results exceed the Ref. Range and clogging still exist after other maintenance procedures have been adopted.
- Analyzer has been running for more than 24 hours.

The cleanser soak procedures are shown as follows.

1. Click the **Service** icon in the menu page to access the **Service** interface.

Figure 12-9 Service



2. Click Maintain in the Maintenance selection.

The interface as shown in the following picture will pop up on the screen.



3. Click the icon of Cleanser Soak.

A dialog box as shown below will pop up.


Figure 12-10 Cleanser Soak



4. Click Yes.

A dialog box as shown below will pop up.


Figure 12-11 Cleanser Soak Prompt



5. Present the cleanser to the sample probe as per the prompt, and press the aspirate key or click the Aspirate button.

"Cleanser soaking..." and the soaking time will appear as shown below. As shown below.

Figure 12-12 Cleanser Soaking Process Prompt




After one minute of soaking, you can stop it manually.

6. Click the **Stop soaking** button, or wait for 19 minutes until the automatic soaking is completed.

After the soaking is completed, a prompt "Cleanser Maintenance done!" will appear. See Figure 12-13.

Figure 12-13 Cleanser Maintenance Done



- 7. Click Close.
- 8. Perform the above procedures to perform the cleanser soak again if necessary.

12.2.4 Comprehensive Device Maintenance


The Comprehensive Device Maintenance feature includes fluidics initialization, comprehensive device cleaning, emptying fluidics and preparing to ship.

12.2.4.1 Fluidics Initialization

After maintaining the fluidic system or replacing a main part of the analyzer, you should perform this procedure to initialize the fluidic system.

Do as follows to perform the fluidics initialization:

1. Click the **Service** icon in the menu page to access the **Service** interface.



2. Click Comprehensive Device in the Maintenance selection.


The interface as shown below will pop up on the screen.

Figure 12-14 Comprehensive Device Maintenance



3. Double click the icon of Fluidics Initialization.

The analyzer starts to perform the fluidics initialization procedure. After the initialization is complete, a message box will pop up.



4. Click OK.

12.2.4.2 Clean Fluidics

If the background results of parameters are out of the background range, the comprehensive device cleaning should be cleansed.

Procedures for comprehensive device cleaning are shown as below:


- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Comprehensive Device in the Maintenance selection.

A dialog box will pop up as shown below.



3. Click the icon of Clean Fluidics.

The analyzer starts to perform the fluidics cleaning procedure. After the cleaning is completed, the following message box will pop up.



4. Click OK.

12.2.4.3 Empty Fluidics

This function enables the device to empty fluidics to prevent crystallization and maintain device performance when the device has not been used for more than one week.

Procedures for emptying fluidics are shown as below:


- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Comprehensive Device in the Maintenance selection.

The interface as shown below will pop up on the screen.



3. Click the icon of Empty Fluidics.

A dialog box will pop up as shown below.



- 4. Click Yes.
- 5. Remove all reagent pickup tube assemblies according to the prompt, and then click **OK** to start emptying the fluidic system.

A dialog box will pop up to prompt you to power off the device after the emptying is complete.

Empty Fluidics done. Please power off the analyzer!

- 6. Turn to [O] the [O/I] switch located on the left side of the analyzer.
- 7. After shutdown, empty the waste in the waste container, and dispose of it.



WARNING

Be sure to dispose of reagents, waste, samples, consumables, etc. according to local legislations and regulations.

12.2.4.4 Prepare to Ship

If the analyzer is not to be used for over two weeks or needs be transported over a long distance (transporting time>2h), you should perform this procedure.

Do as follows to perform the prepare-to-ship procedure:


- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Comprehensive Device in the Maintenance selection.

The interface as shown below will pop up on the screen.




3. Click the icon of Prepare to Ship.

A message box shown below will pop up.




4. Click Yes.

A dialog box will pop up as shown below.



5. Remove all reagent pickup tube assemblies according to the prompt, and then click **OK** to start emptying the fluidic system.

After the emptying is complete, a message box will pop up.




6. Place all reagent pickup tube assemblies into the distilled water, and then click **OK** to start priming.

NOTE

- Be sure to use distilled water in order to ensure the normal use of the device in the future. In addition, the beaker holding the distilled water needs to be cleaned thoroughly.
- The diluent pipe and lyse pipes should be stored separately in two beakers.

System performs the filling operation. After the filling is completed, the following dialog box will pop up.



7. Take out the diluent and lyse pipes from the distilled water as per the prompt, then click **OK**. After the operation, dialog box will pop up to prompt you to power off the device.

Prepare to Ship done. Please power off the analyzer!

- 8. Turn to [O] the [O/I] switch located on the left side of the analyzer.
- 9. After shutdown, empty the waste in the waste container, and dispose of it.



WARNING

Be sure to dispose of reagents, waste, samples, consumables, etc. according to local legislations and regulations.

12.2.5 Auto Clean

There will be a certain amount of contamination accumulated after running a certain amount of samples without shutting down the analyzer. When the sample count amounts to over 100, the analyzer will perform the cleaning procedure automatically once, and a prompt will be displayed on the screen.

In addition, the analyzer will perform the auto clean procedures if there has been no fluidics sequential operation for more than one hour.

NOTE

Once the auto clean is performed or the analyzer is shut down, the statistical data of auto clean will be cleared automatically.

12.2.6 Auto Prompt for Cleanser Soak

If the analyzer has been running for more than 24 hours but hasn't performed cleanser maintenance when the auto maintenance time is reached, the system will prompt to perform cleanser soak immediately, so as to prevent the accumulation of contamination.

Click **Yes**, then you can perform the cleanser maintenance as per the prompt and the description in **12.2.3.2 Cleanser Soak**.

Click **No**, then the system will remind you every 10 minutes until you perform the maintenance.

NOTE

- At the Self-test or Status interface, the analyzer does not ask for confirmation to perform the cleanser soak.
- If the analyzer is running or has problems when the conditions of auto prompt for cleanser soak
 is satisfied, the analyzer will prompt again after the current operation is completed or the
 problems are resolved.
- After cleanser soak is completed, the accumulative count values will be cleared automatically.
- Cleanser soak is an important step in comprehensive device maintenance. It is recommended not to stop soaking halfway.

12.2.7 Auto Sleep

When the fluidics system stops working for 60 minutes (default setting), then the analyzer will enter the sleeping status automatically. You can change the waiting time for auto sleeping as needed, see 5.3.4 Auto Maintenance.

When the analyzer is in the sleep mode, a prompt will be displayed on the screen. Touch the screen or press the aspirate key on the analyzer to wake it up.

NOTE

- If it is the time to auto sleep but the analyzer is error status, then only after the error is removed will auto sleep start accordingly.
- Different maintenances will be performed by the analyzer automatically when exiting the sleep mode, and the exiting time depends on how long the analyzer was in the sleep mode.
- If errors occur when you are trying to cancel the auto sleep of the analyzer, please refer to 13
 Troubleshooting for solving the problems.

12.3 Self-inspection

This feature is to test if some important components of the device can function properly or not, including syringe and sampling assembly self-inspection, pressure and vacuum self-inspection, valve self-inspection and other self-inspections.

NOTE

If the testing result is abnormal, you should try again for several times; if the abnormalities persist, please contact MR customer service department or your local agent.

12.3.1 Syringe and Sampling Mechanism

You can test the performance of all syringes and sampling mechanisms.

The self-inspection procedures are shown as below:

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Syringe Self-test in the Self-test selection.

The interface as shown in Figure 12-15 will pop up on the screen.


Figure 12-15 Syringe



3. Double click the part that needs to be tested, e.g. **Sample Syringe**, and wait for the self-inspection results.

After the self-test is completed, a dialog box will pop up to show the self-test results.

Figure 12-16 Syringe Self-test Results



4. Click **OK** to close the message box.

12.3.2 Pressure and Vacuum


This feature is to test the pressure and vacuum inside the device.

Procedures for pressure (or vacuum) self-inspection are shown as below:

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Pressure Self-test in the Self-test selection.


The interface as shown in Figure 12-17 will pop up on the screen.

Figure 12-17 Pressure and Vacuum Self-inspection



3. Click Pressure (or Vacuum).

The system will perform the corresponding self-test operations. After the self-test is completed, a dialog box will pop up to show the self-test results.



4. Click **OK** to close the message box.

12.3.3 Valve & Pump

When controlling the switches of different valves (pumps), you can judge if the valves (pumps) are operating properly by the sound of opening, closing or manually touching the corresponding valves (pumps).

The procedures for valve self-inspection are shown as follows:

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Valve/Pump Self-test in the Self-test selection.

The interface as shown in Figure 12-18 will pop up on the screen.

Valve:

1 2 3 4 5 6

7 8 9 10 11

Pump:

Figure 12-18 Valve/Pump Self-test

3. Click the desired Valve No. (e.g. 1), then confirm whether it works properly by the sound of its opening and closing.

12.3.4 Others

You can perform the self-test for WBC and RBC aperture voltage.

Take the self-test for RBC aperture voltage as an example, its operation steps are as follows:

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Other Self-test in the Self-test selection.

The interface as shown in Figure 12-19 will pop up on the screen.


Figure 12-19 Other Self-test



3. Click RBC Aperture Voltage to start self-test.

The system will perform the corresponding self-test operations. After the self-inspection is completed, a dialog box will pop up to show the self-inspection results.

Figure 12-20 RBC Aperture Voltage Self-test Results



12.4 System Status


You can view the current status information of the analyzer in the **Status** selection, including temperature, voltage and current, and version information.

12.4.1 Temperature

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click **Temperature** in the **Status** selection.

The interface as shown in Figure 12-21 will pop up on the screen.

Figure 12-21 View Temperature Status




You can view the current ambient temperature of the analyzer. If the results of the temperature testing exceed the normal range, they will be highlighted by the red background.

12.4.2 Voltage and Current

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Voltage/Current.

The interface as shown below will pop up on the screen.

Figure 12-22 Voltage and Current



You can view the voltage and current information of the analyzer. The voltage or current value that exceeds the normal range will be displayed in a red background.


12.4.3 Version Information

You can view the current version information of all parts of the analyzer, and export the version information to a USB flash disk. Specific steps are shown below:

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click **Version Info.** In the **Status** selection.


Version information interface will pop up on the screen. See Figure 12-23.

Figure 12-23 Version Information




- 3. Insert a USB flash disk in the USB interface on the analyzer.
- 4. Click **Export**, and select the export path in the dialog box, and then enter the file name.

The file will be exported to the root directory of the USB flash disk (/udisk/sda1) by default as shown below.



5. Click **Save** to start exporting.

After Export is completed, the message box as shown below will pop up.



6. Click OK to exit.

12.5 Log

In the Log selection, you can view the records of All Logs, Set Paras, Fault Logs and Other Logs.

NOTE

- If a new record is added when the log is full, the newest record will overwrite the oldest one automatically.
- The administrator can view both his/her own operation logs and the general users' operation logs, while the general users can only review their own operation logs.
- The log can keep records of up to 5 years.

12.5.1 All Logs

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click **All Logs** in the **Log** selection.

You can view all logs (visible to the users of the current access level).




Figure 12-24 All Logs

3. Select the dates in the two date textboxes, and then you can view the all logs within the date range, including operation time, log information and the operator.


12.5.2 Parameter Revision Logs

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Set Paras in the Log selection.

You can view the parameter revision logs (which can be viewed by the user with the current level

of access) within a specified date range.

Figure 12-25 Parameter Revision Logs



3. Select the dates in the two date textboxes, and then you can view the parameter revision logs within the date range, including the revision date and time, revision summary and the operator.

12.5.3 Fault Logs

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Fault Logs in the Log selection.

You can view all logs (visible to the users of the current access level).




Figure 12-26 Fault Logs

3. Select the dates in the two date textboxes, and then you can view the fault logs within the date range, including date and time when the faults occur, fault description and the operator.

12.5.4 Other Logs

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Other Logs in the Log selection.

You can view other logs besides parameter revision logs and fault logs.

Other Logs 2015/07/23 2015/07/23 **Summary Information** Details Operator No. Time Background mode counting run succ 2015/07/23 18:55:21 Run 2 Background mode counting run succ... Administrator a... 3 2015/07/23 18:55:17 Startup Startup Administrator a... 2015/07/23 18:55:15 Login 4 admin(admin) Login Administrator a... 5 2015/07/23 18:06:38 Shutdown after Pre... Shutdown after Prepare to Ship proce.. Administrator a... 6 2015/07/23 18:06:38 Run Background mode counting run succ... Administrator a... 7 2015/07/23 18:06:35 Run Background mode counting run succ... Administrator a... 8 2015/07/23 18:06:31 Startup Startup Administrator a... 9 2015/07/23 18:06:29 Login admin(admin) Login Administrator a.. 10 2015/07/23 18:04:42 | Shutdown after Pre... Shutdown after Prepare to Ship proce.. Administrator a... Date and Time: 2015/07/23 18:55:24 Operator: Administrator admin (admin) Summary Information: Run Details: Background mode counting run successfully

Figure 12-27 Other Logs

3. Select the dates in the two date textboxes to view the logs within the date range, including operation date and time, operation records and the operator.

12.6 Touch Screen Calibration

When the touch screen has offset, it needs to be recalibrated. Specific steps are shown below:

- 1. Click the **Service** icon in the menu page to access the **Service** interface.
- 2. Click Touch Screen Cal. in the Cal selection.
- 3. Click the calibration point "+" on the screen in order.

When the calibration point disappears and the system return to the service screen, it indicates the completion of the calibration.

13 Troubleshooting

13.1 Introduction

This chapter contains information that is helpful in locating and resolving problems that may occur during the operation of your analyzer.

NOTE

This chapter is not a complete service manual and is limited to problems that are readily diagnosed and/or corrected by the user of the analyzer. If the recommended solution fails to solve the problem, contact MR customer service department or your local agent.

13.2 Dealing with Error Messages

In the use of the analyzer, when the software detects abnormalities, an error message will be displayed on the upper right of the screen as shown in Figure 13-1 and the main unit will sound an alarm.

Figure 13-1 Error Messages

Background abnormal.

You can refer to the following steps to deal with the error messages.

1. Click the error message area.

As shown in Figure 13-2, the popup dialog box displays the error description and its help information. The error descriptions are displayed in the order of error occurrence.




Figure 13-2 Error Message Dialog Box

- 2. Touch the screen to disable the beep.
- 3. Click Remove Error.

Normally, the system will automatically remove the errors and close the dialog box.

For errors which cannot be removed automatically, you can take appropriate actions by following the error help information or **13.3** *Error Message Reference*.

13.3 Error Message Reference

Possible errors and the corresponding help information are shown in Table 13-1.

Table 13-1 Error Message Reference

Error description	Troubleshooting Information	
Abnormal -12V power.	Please power off the analyzer directly and restart later. If the error still exists, contact our customer service department.	
Abnormal voltage of constant-current voltage abnormal.	Please power off the analyzer directly and restart later. If the error still exists, contact our customer service department.	
Startup failure.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Startup initialization is not executed.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	

Error description	Troubleshooting Information	
Right side door is open.	 Close the right side door. Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
Abnormal +12V power.	Please power off the analyzer directly and restart later. If the error still exists, contact our customer service department.	
	Check if the DIL-E diluent expires. If so, replace it with a new container of DIL-E.	
DIL-E expiration.	2. Click the Remove Error button, the Reagent Management screen will be displayed.	
	3. Set the reagent information by referring to <i>11 Reagent Management</i> .	
	4. If the error still exists, contact our customer service department.	
	1. Check if the LYE-1 lyse expires. If so, replace it with a new container of LYE-1.	
LYE-1 expiration.	2. Click the Remove Error button, the Reagent Management screen will be displayed.	
	3. Set the reagent information by referring to <i>11 Reagent Management</i> .	
	4. If the error still exists, contact our customer service department.	
Abnormal HGB background voltage	1. Adjust the HGB background voltage within the specified range (4.2V~4.8V), preferably 4.5V. Refer to <i>5.5.1 Gain Settings</i> .	
voltage	2. If the error still exists, contact our customer service department.	
Abnormal RBC aperture voltage.	 Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
Abnormal WBC aperture voltage	 Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
	Click the Remove Error button to remove this error.	
RBC clogging	2. If the error is reported frequently, see <i>12.2.3.2 Cleanser Soak</i> to dip the RBC bath in the cleanser.	
	3. If the error still exists, contact our customer service department.	
	Click the Remove Error button to remove this error.	
WBC clogging.	2. If the error is reported frequently, see <i>12.2.3.2 Cleanser Soak</i> to dip the RBC bath in the cleanser.	
	3. If the error still exists, contact our customer service department.	
	Check whether the diluent is contaminated.	
Abnormal background	If not, click the Remove Error button to remove the error.	
	If the error still exists, contact our customer service department.	
Failed to read sample	Click the Remove Error button to remove this error.	
syringe parameter.	2. If the error still exists, contact our customer service department.	

Error description	Troubleshooting Information	
Failed to configure sample syringe parameter.	 Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
Sample syringe timeout	 Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
Sample syringe is busy.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Command parameter error of the sampling assembly.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Sampling assembly timeout	 Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
Sampling assembly is busy.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Vertical motor instruction parameter error.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Failed to read vertical motor parameter.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Vertical motor timeout	 Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
Failed to read the remaining steps of vertical motor.	1. Click the Remove Error button to remove this error. 2. If the error still exists, contact our customer service department.	
Vertical motor is busy.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Failed to read ambient temperature.	Make sure the temperature sensor is correctly installed. If the error still exists, contact our customer service department.	
1. Empty the waste container or install a new waste container or unstall a new waste container or install a new waste container or i		
Failed to read horizontal motor parameter.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Failed to configure Horizontal motor parameter.	1. Click the Remove Error button to remove this error. 2. If the error still exists, contact our customer service department.	
Horizontal motor timeout	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	
Abnormal horizontal-motor photocoupler.	Click the Remove Error button to remove this error. If the error still exists, contact our customer service department.	

Error description	Troubleshooting Information	
Horizontal motor is busy.	 Click the Remove Error button to remove this error. If the error still exists, contact our customer service department. 	
	Check whether the DIL-E container is empty. If so, perform step 2; or if there is still plenty of DIL-E contact our customer service department. Install a new container of DIL-E.	
No DIL-E.	3. Click the Remove Error button to remove this error.	
	If the error still exists after a new container of DIL-E is installed, contact our customer service department.	
No LYE-1.	Check whether the LYE-1 container is empty. If so, perform step 2; or if there is still plenty of LYE-1, contact our customer service department. Install a new container of LYE-1.	
	3. Click the Remove Error button to remove this error.	
	4. If the error still exists after a new container of LYE-1 is installed, contact our customer service department.	
DIL-E not replaced.	Click the Remove Error button to remove this error.	
DIL-E not replaced.	2. If the error still exists, contact our customer service department.	
LYE-1 not replaced.	Click the Remove Error button to remove this error.	
ETE THOUTOplaced.	2. If the error still exists, contact our customer service department.	
Abnormal 12V driving	Please power off the analyzer directly and restart later.	
power supply.	2. If the error still exists, contact our customer service department.	
Abnormal 24V driving	Please power off the analyzer directly and restart later.	
power supply.	2. If the error still exists, contact our customer service department.	
	Check whether the LYE-1 container is empty. If so, perform step 2; or if there is still plenty of LYE-1, contact our customer service department.	
	2. Install a new container of LYE-1.	
Insufficient LYE-1.	3. Click the Remove Error button, the Reagent Management screen will be displayed.	
	4. Set the reagent information by referring to <i>11 Reagent Management</i> .	
	5. If the error still exists, contact our customer service department.	
	1. Check whether the DIL-E container is empty. If so, perform step 2; or if there is still plenty of DIL-E, contact our customer service department.	
	2. Install a new container of DIL-E.	
Insufficient DIL-E.	3. Click the Remove Error button, the Reagent Management screen will be displayed.	
	4. Set the reagent information by referring to <i>11 Reagent Management</i> .	
	5. If the error still exists, contact our customer service department.	

Appendix A Specifications

A.1 Classification

According to the CE classification, the Auto Hematology Analyzer belongs to in vitro diagnostic medical devices, rather than those covered by Annex II and devices for performance evaluation.

A.2 Reagents

Reagent Type	Reagent Name
Diluent	DIL-E Diluent
Lyse	LYE-1 LYSE
Medical cleanser	CLE-P Cleanser

A.3 Parameters

Parameter	Abbreviation	Default Unit
White Blood Cell count	WBC	10 ⁹ /L
Number of Granulocytes	Gran#	10 ⁹ /L
Number of Lymphocytes	Lym#	10 ⁹ /L
Number of Mid-sized Cells	Mid#	10 ⁹ /L
Percentage of Granulocytes	Gran%	%
Percentage of Lymphocytes	Lym%	%
Percentage of Mid-sized Cells	Mid%	%
Red Blood Cell count	RBC	10 ¹² /L
Hemoglobin Concentration	HGB	g/L
Hematocrit	НСТ	%
Mean Corpuscular Volume	MCV	fL
Mean Corpuscular Hemoglobin	MCH	pg
Mean Corpuscular Hemoglobin Concentration	MCHC	g/L

Parameter	Abbreviation	Default Unit
Red Blood Cell Distribution Width Standard Deviation (RDW-SD)	RDW-SD	fL
Red Blood Cell Distribution Width Coefficient of Variation (RDW-CV)	RDW-CV	%
Platelet count (PLT count, 10 ⁹ /L)	PLT	10 ⁹ /L
Mean Platelet Volume (MPV, fL)	MPV	fL
Platelet Distribution Width (PDW)	PDW	N/A
Plateletcrit (PCT)	PCT	%
Platelet-large cell ratio	P-LCR	%
Platelet-large cell count	P-LCC	10 ⁹ /L
Red Blood Cell Histogram	RBC Histogram	N/A
Platelet Histogram	PLT Histogram	N/A
White Blood Cell Histogram	WBC Histogram	N/A

A.4 Performance Specifications

A.4.1 Display Range

Parameter	Linearity Range	Display Range
WBC	0~300×10 ⁹ /L	0~999×10 ⁹ /L
RBC	0.00~8.50×10 ¹² /L	0~18.00×10 ¹² /L
HGB	0~250g/L	0~300g/L
PLT	0~3000×10 ⁹ /L	0~5000×10 ⁹ /L
НСТ	0~67%	0%~80%

A.4.2 Normal Background

Parameter	Normal Background
WBC	≤0.2×10 ⁹ /L
RBC	≤0.02×10 ¹² /L
HGB	≤1g/L
PLT	≤5×10 ⁹ /L
нст	≤0.5%

A.4.3 Linearity Range

Parameter	Linearity range	Deviation range (Whole blood mode)
WDC	(0.00~100.00)×10 ⁹ /L	±0.30×109/L or ±5%
WBC	(100.01~300.00)×10 ⁹ /L	±10%
RBC	(0.00~8.50)×10 ¹² /L	±0.05×1012/L or ±5%
HGB	(0~250) g/L	±2g/L or ±2%
DLT	(0~1000)×10 ⁹ /L (RBC≤7.0)	±1 0×109/L or ±8%
PLT 1001~3000×10 ⁹ /L (RBC≤7.0		±12%
HCT	0~67%	±2% (HCT value) or ±3% (deviation percent)

A.4.4 Repeatability

These repeatability requirements apply only to the situation in which a qualified sample has been run for 11 times and the results of the 2nd to 11th runs are used to calculate the repeatabilities.

Parameter	Condition	Repeatability (CV%/absolute deviation d*)
WBC	(7.0~15.0)×10 ⁹ /L	≤2.0%
	(4.0~6.9)×10 ⁹ /L	≤3.5%
RBC	(3.5~6.5)×10 ¹² /L	≤1.5%
HGB	(100~180)g/L	≤1.5%
MCV	(70~110)fL	≤0.5%
PLT	(150~500)×10 ⁹ /L	≤4.0%

Absolute deviation d = analysis result – average of analysis results

A.4.5 Carryover

Parameter	Carryover
WBC	≤0.5%
RBC	≤0.5%
HGB	≤0.5%
PLT	≤1.0%
WBC	≤0.5%

A.5 Input/output Device



WARNING

Accessory equipment connected to the analogue and digital interfaces must comply with the relevant Safety and EMC standards (e.g., IEC 60950 Safety of Information Technology Equipment Standard and CISPR 22 EMC of Information Technology Equipment Standard (CLASS B)). Anyone who connects additional equipment to the signal input or output ports and configures an IVD system is responsible for ensuring that the system works properly and complies with the safety and EMC requirements. If you have any problem, consult the technical services department of your local agent.

Host

- Touch screen: 10.4 inches embedded touch screen with a resolution of 800×600
- > Thermal printer
- One LAN interface
- 4 USB interfaces
- Power
 - ➤ Voltage: A.C 100V~240V
 - ➤ Input power: ≤200VA
 - > Frequency: 50/60 Hz
- Keyboard (Optional, USB)
- Mouse (Optional, USB)
- External barcode scanner (optional, USB)
- Printer (optional, USB)
- USB flash disk (optional, USB)

A.6 EMC Description

This equipment complies with the emission and immunity requirements of the IEC 61326-1:2012, EN 61326-1:2013. IEC 61326-6-2-6:2012 and EN 61326-2-6:2013.

This equipment has been designed and tested to CISPR 11 Class A. In a domestic environment it may cause radio interference, in which case, you may need to take measures to mitigate the interference.


A.7 Environmental conditions



Be sure to use and store the analyzer in the specified environment.

Environmental conditions	Operating Environment	Storage Environment	Running Environment
Ambient temperature	10°C~30°C	-10°C~40°C	5°C~40°C
Relative humidity	20%~85%	10%~90%	10%~90%
Atmospheric pressure	70kPa~106kPa	50kPa~106kPa	70kPa~106kPa

A.8 Dimensions and Weight



Analyzer	Dimensions and Weight
Width (mm)	≈360
Height (mm)	≈455
Depth (mm)	≈410
Weight (kg)	≈24

A.9 Expected service life

More than 5 years.

A.10 Contraindications

N/A

Appendix B Terms and Abbreviations

CWB Capillary Whole Blood

PD Predilute Blood RF

Radio Frequency

VWB Venous Whole Blood

Appendix C Packing List

No.	Name	Quantity
1	Auto Hematology Analyzer	1
2	Power cable	1
3	Grounding cable	1
4	Operator's Manual	1
5	Software installation CD-ROM	1
6	Quick Operation Guide Card	1
7	Diluent Adapter Tube	1
8	Waste Float Adapter Tube	1
9	DIL-E Diluent (20L)	1
10	LYE-1 Lyse (200mL)	1
11	Cleanser (50mL)	1
12	Warranty Card	1
13	Waste container	1
14	Auto Hematology Analyzer Inspection Record	1
15	Certificate of Delivery and Acceptance	1
16	Packing list	1
17	Reagent Operation Guide for Closed System (for closed system only)	1

P/N: 65020020A (2.Q)